

UniDec

Generated by Doxygen 1.8.10

Tue Sep 22 2015 16:47:44

Contents

1	Hierarchical Index	1
1.1	Class Hierarchy	1
2	Class Index	3
2.1	Class List	3
3	Class Documentation	5
3.1	miscwindows.AdditionalParameters Class Reference	5
3.2	PlotAnimations.AnimationWindow Class Reference	6
3.3	masstools.AutocorrWindow Class Reference	7
3.4	plot1d.BarChart Class Reference	7
3.5	nativez.ColorList Class Reference	8
3.5.1	Detailed Description	8
3.5.2	Constructor & Destructor Documentation	9
3.5.2.1	__init__(self, parent)	9
3.6	ColorPlot.ColorPlot2D Class Reference	9
3.7	masstools.CorrListCtrl Class Reference	10
3.8	masstools.CorrListCtrlPanel Class Reference	10
3.9	plot3d.CubePlot Class Reference	10
3.10	datacollector.DataCollector Class Reference	12
3.11	unidec.DataContainer Class Reference	14
3.11.1	Constructor & Destructor Documentation	15
3.11.1.1	__init__(self)	15
3.12	Extract2D.Extract2DPlot Class Reference	15
3.13	miscwindows.FileNameDialog Class Reference	16
3.14	import_wizard.ImportWizard Class Reference	18
3.14.1	Member Function Documentation	18
3.14.1.1	export_file(self, evt)	18
3.14.1.2	export_then_load(self, evt)	18
3.14.1.3	get_folder_path(self, evt)	18
3.15	IM_windows.IMToolExtract Class Reference	19
3.16	IM_windows.IMTools Class Reference	20

3.17	UniFit.KDmodel Class Reference	21
3.18	UniFit.kdstruct Class Reference	22
3.19	datacollector.ListCtrlPanel Class Reference	23
3.20	mainwindow.Mainwindow Class Reference	23
3.20.1	Detailed Description	27
3.20.2	Constructor & Destructor Documentation	27
3.20.2.1	__init__(self, parent, title, config)	27
3.20.3	Member Function Documentation	27
3.20.3.1	clear_all_plots	27
3.20.3.2	export_gui_to_config(self)	27
3.20.3.3	import_config_to_gui(self)	28
3.20.3.4	menu_401_403(self, event)	28
3.20.3.5	on_about(self, e)	28
3.20.3.6	on_check_manual(self, e)	28
3.20.3.7	on_defaults(self, e)	28
3.20.3.8	on_exit(self, e)	28
3.20.3.9	on_flip_mode(self, e)	28
3.20.3.10	on_flip_tabbed(self, e)	28
3.20.3.11	on_flip_twave(self, e)	28
3.20.3.12	on_mass_list(self, e)	29
3.20.3.13	on_motion(self, xpos, ypos)	29
3.20.3.14	on_save_figure_dialog(self, e)	29
3.20.3.15	on_save_figure_eps(self, e)	29
3.20.3.16	on_save_figure_pdf(self, e)	29
3.20.3.17	on_save_figure_png(self, e, kwargs)	29
3.20.3.18	on_save_figure_small(self, e)	29
3.20.3.19	save_all_figures(self, extension, extension2="", e=0, header=None, kwargs)	29
3.20.3.20	setup_main_panel(self)	30
3.20.3.21	setup_menu(self)	30
3.20.3.22	setup_shortcuts(self)	30
3.20.3.23	setup_tool_tips(self)	30
3.20.3.24	shrink_all_figures(self)	30
3.20.3.25	shrink_figure(self, plot)	30
3.21	masstools.ManualSelection Class Reference	30
3.22	MassModelFitter.mass Class Reference	31
3.23	MassDefects.MassDefectWindow Class Reference	32
3.24	MassFitter.MassFitter Class Reference	33
3.25	masstools.MassSelection Class Reference	34
3.26	mzMLImporter.mzMLImporter Class Reference	36
3.27	nativez.NativeZ Class Reference	36

3.28	datacollector.NetworkFrame Class Reference	38
3.29	peakstructure.Peak Class Reference	38
3.29.1	Detailed Description	39
3.29.2	Constructor & Destructor Documentation	39
3.29.2.1	__init__(self)	39
3.30	peakstructure.Peaks Class Reference	39
3.30.1	Detailed Description	40
3.30.2	Constructor & Destructor Documentation	40
3.30.2.1	__init__(self)	40
3.30.3	Member Function Documentation	40
3.30.3.1	add_peaks(self, parray)	40
3.30.3.2	default_params	40
3.30.3.3	get_mass_defects	40
3.30.3.4	score_peaks	40
3.31	peakwidthtools.PeakTools1d Class Reference	40
3.32	peakwidthtools.PeakTools2d Class Reference	42
3.33	twitter_interface.PinWindow Class Reference	43
3.34	plot1d.Plot1d Class Reference	44
3.35	plot2d.Plot2d Class Reference	44
3.36	plot2d.Plot2dMass Class Reference	45
3.37	PlottingWindow.PlottingWindow Class Reference	46
3.38	miscwindows.SaveFigureDialog Class Reference	47
3.39	GUniDec.Shell Class Reference	48
3.40	miscwindows.SingleInputDialog Class Reference	48
3.41	peaklistsort.TestListCtrl Class Reference	49
3.42	masstools.TestListCtrl Class Reference	49
3.43	masstools.TestListCtrl2 Class Reference	49
3.44	masstools.TestListCtrl3 Class Reference	50
3.45	IM_windows.TestListCtrl4 Class Reference	51
3.46	masstools.TestListCtrlMatch Class Reference	51
3.47	masstools.TestListCtrlPanel Class Reference	52
3.48	peaklistsort.TestListCtrlPanel Class Reference	52
3.49	masstools.TestListCtrlPanel2 Class Reference	53
3.50	masstools.TestListCtrlPanel3 Class Reference	53
3.51	IM_windows.TestListCtrlPanel4 Class Reference	54
3.52	masstools.TestListCtrlPanelMatch Class Reference	55
3.53	import_wizard_treectrl.TreeCtrl Class Reference	55
3.54	import_wizard_treectrl.TreeCtrlPanel Class Reference	55
3.54.1	Member Function Documentation	56
3.54.1.1	add_children	56

3.54.1.2	add_root	56
3.55	twitter_interface.TwitterWindow Class Reference	56
3.56	unidec.UniDec Class Reference	58
3.56.1	Constructor & Destructor Documentation	59
3.56.1.1	__init__(self)	59
3.56.2	Member Function Documentation	59
3.56.2.1	autocorrelation	59
3.56.2.2	autointegrate	59
3.56.2.3	center_of_mass	59
3.56.2.4	check_badness(self)	59
3.56.2.5	convolve_peaks(self)	60
3.56.2.6	cross_validate	60
3.56.2.7	export_config	60
3.56.2.8	export_params(self, e)	60
3.56.2.9	fit_all_masses(self)	60
3.56.2.10	get_charge_peaks(self)	60
3.56.2.11	initialize(self)	60
3.56.2.12	integrate	60
3.56.2.13	kendrick_continuous	61
3.56.2.14	kendrick_peaks	61
3.56.2.15	load_config(self, f_name)	61
3.56.2.16	load_default(self)	61
3.56.2.17	load_state(self, load_path)	61
3.56.2.18	mass_grid_to_f_grid(self)	61
3.56.2.19	normalize_peaks(self)	61
3.56.2.20	open_file(self, file_name, file_directory, kwargs)	62
3.56.2.21	pick_peaks(self)	62
3.56.2.22	process_data(self, kwargs)	62
3.56.2.23	process_mass_data(self)	62
3.56.2.24	raw_process	62
3.56.2.25	reset_config(self)	62
3.56.2.26	run_unidec	63
3.56.2.27	save_default(self)	63
3.56.2.28	unidec_imports	63
3.57	GUniDec.UniDecApp Class Reference	63
3.57.1	Detailed Description	65
3.57.2	Member Function Documentation	65
3.57.2.1	export_config	65
3.57.2.2	import_config	65
3.58	unidecstructure.UniDecConfig Class Reference	65

3.58.1 Detailed Description	67
3.58.2 Constructor & Destructor Documentation	68
3.58.2.1 __init__(self)	68
3.58.3 Member Function Documentation	68
3.58.3.1 check_badness(self)	68
3.58.3.2 config_export(self, name)	68
3.58.3.3 config_import(self, name)	68
3.58.3.4 default_colormaps(self)	68
3.58.3.5 default_file_names(self)	68
3.58.3.6 default_high_res(self)	68
3.58.3.7 default_zero_charge(self)	68
3.58.3.8 initialize(self)	69
3.58.3.9 initialize_system_paths(self)	69
3.58.3.10 print_config(self)	69
3.59 unidec_tests.UniDecTest Class Reference	69
3.60 import_wizard_grid.WizardGrid Class Reference	70
3.60.1 Detailed Description	70
3.60.2 Member Function Documentation	70
3.60.2.1 fill_down(self, evt)	70
3.60.2.2 set_labels(self, mode)	70
3.60.2.3 showPopupMenu(self, evt)	71
3.61 datacollector.XListCtrl Class Reference	71
3.62 datacollector.YListCtrl Class Reference	71
3.63 nativez.zoffset Class Reference	72
3.64 ZoomBox.ZoomBox Class Reference	72
3.64.1 Detailed Description	73
3.64.2 Constructor & Destructor Documentation	74
3.64.2.1 __init__	74
3.64.3 Member Function Documentation	74
3.64.3.1 get_active(self)	74
3.64.3.2 set_active(self, active)	74
3.65 ZoomSpan.ZoomSpan Class Reference	74
3.65.1 Detailed Description	75
3.65.2 Constructor & Destructor Documentation	76
3.65.2.1 __init__	76
Index	77

Chapter 1

Hierarchical Index

1.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

ColumnSorterMixin	
peaklistsort.TestListCtrlPanel	52
unidec.DataContainer	14
Dialog	
IM_windows.IMToolExtract	19
IM_windows.IMTools	20
masstools.AutocorrWindow	7
masstools.ManualSelection	30
masstools.MassSelection	34
miscwindows.AdditionalParameters	5
miscwindows.FileNameDialog	16
miscwindows.SaveFigureDialog	47
miscwindows.SingleInputDialog	48
nativez.NativeZ	36
peakwidthtools.PeakTools1d	40
peakwidthtools.PeakTools2d	42
twitter_interface.PinWindow	43
twitter_interface.TwitterWindow	56
Frame	
datacollector.DataCollector	12
Extract2D.Extract2DPlot	15
import_wizard.ImportWizard	18
mainwindow.Mainwindow	23
MassDefects.MassDefectWindow	32
PlotAnimations.AnimationWindow	6
Grid	
import_wizard_grid.WizardGrid	70
UniFit.KDmodel	21
UniFit.kdstruct	22
ListCtrl	
datacollector.XListCtrl	71
datacollector.YListCtrl	71
IM_windows.TestListCtrl4	51
masstools.CorrListCtrl	10
masstools.TestListCtrl	49
masstools.TestListCtrl2	49
masstools.TestListCtrl3	50
masstools.TestListCtrlMatch	51

peaklistsort.TestListCtrl	49
ListCtrlAutoWidthMixin	
datacollector.XListCtrl	71
datacollector.YListCtrl	71
IM_windows.TestListCtrl4	51
masstools.CorrListCtrl	10
masstools.TestListCtrl	49
masstools.TestListCtrl2	49
masstools.TestListCtrl3	50
masstools.TestListCtrlMatch	51
MassModelFitter.mass	31
MassFitter.MassFitter	33
mzMLImporter.mzMLImporter	36
object	
GUniDec.Shell	48
GUniDec.UniDecApp	63
Panel	
datacollector.ListCtrlPanel	23
IM_windows.TestListCtrlPanel4	54
import_wizard_treectrl.TreeCtrlPanel	55
masstools.CorrListCtrlPanel	10
masstools.TestListCtrlPanel	52
masstools.TestListCtrlPanel2	53
masstools.TestListCtrlPanel3	53
masstools.TestListCtrlPanelMatch	55
nativez.ColorList	8
peaklistsort.TestListCtrlPanel	52
peakstructure.Peak	38
peakstructure.Peaks	39
TextEditMixin	
datacollector.XListCtrl	71
datacollector.YListCtrl	71
IM_windows.TestListCtrl4	51
masstools.CorrListCtrl	10
masstools.TestListCtrl	49
masstools.TestListCtrl2	49
masstools.TestListCtrl3	50
masstools.TestListCtrlMatch	51
TreeCtrl	
import_wizard_treectrl.TreeCtrl	55
unidec.UniDec	58
unidec_tests.UniDecTest	69
unidecstructure.UniDecConfig	65
Window	
PlottingWindow.PlottingWindow	46
ColorPlot.ColorPlot2D	9
datacollector.NetworkFrame	38
plot1d.BarChart	7
plot1d.Plot1d	44
plot2d.Plot2d	44
plot2d.Plot2dMass	45
plot3d.CubePlot	10
nativez.zoffset	72
ZoomBox.ZoomBox	72
ZoomSpan.ZoomSpan	74

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

miscwindows.AdditionalParameters	5
PlotAnimations.AnimationWindow	6
masstools.AutocorrWindow	7
plot1d.BarChart	7
nativez.ColorList	8
ColorPlot.ColorPlot2D	9
masstools.CorrListCtrl	10
masstools.CorrListCtrlPanel	10
plot3d.CubePlot	10
datacollector.DataCollector	12
unidec.DataContainer	14
Extract2D.Extract2DPlot	15
miscwindows.FileNameDialog	16
import_wizard.ImportWizard	18
IM_windows.IMToolExtract	19
IM_windows.IMTools	20
UniFit.KDmodel	21
UniFit.kdstruct	22
datacollector.ListCtrlPanel	23
mainwindow.Mainwindow	23
masstools.ManualSelection	30
MassModelFitter.mass	31
MassDefects.MassDefectWindow	32
MassFitter.MassFitter	33
masstools.MassSelection	34
mzMLImporter.mzMLImporter	36
nativez.NativeZ	36
datacollector.NetworkFrame	38
peakstructure.Peak	38
peakstructure.Peaks	39
peakwidthtools.PeakTools1d	40
peakwidthtools.PeakTools2d	42
twitter_interface.PinWindow	43
plot1d.Plot1d	44
plot2d.Plot2d	44
plot2d.Plot2dMass	45
PlottingWindow.PlottingWindow	46
miscwindows.SaveFigureDialog	47


GUniDec.Shell	48
miscwindows.SingleInputDialog	48
peaklistsort.TestListCtrl	49
masstools.TestListCtrl	49
masstools.TestListCtrl2	49
masstools.TestListCtrl3	50
IM_windows.TestListCtrl4	51
masstools.TestListCtrlMatch	51
masstools.TestListCtrlPanel	52
peaklistsort.TestListCtrlPanel	52
masstools.TestListCtrlPanel2	53
masstools.TestListCtrlPanel3	53
IM_windows.TestListCtrlPanel4	54
masstools.TestListCtrlPanelMatch	55
import_wizard_treectrl.TreeCtrl	55
import_wizard_treectrl.TreeCtrlPanel	55
twitter_interface.TwitterWindow	56
unidec.UniDec	58
GUniDec.UniDecApp	63
unidecstructure.UniDecConfig	65
unidec_tests.UniDecTest	69
import_wizard_grid.WizardGrid	70
datacollector.XListCtrl	71
datacollector.YListCtrl	71
nativez.zoffset	72
ZoomBox.ZoomBox	72
ZoomSpan.ZoomSpan	74

Chapter 3

Class Documentation

3.1 miscwindows.AdditionalParameters Class Reference

Inheritance diagram for miscwindows.AdditionalParameters:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, config)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)

Public Attributes


- **config**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **hbox5**
- **inputbox5**
- **hbox6**
- **inputbox6**
- **hbox6b**
- **inputbox6b**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/miscwindows.py

3.2 PlotAnimations.AnimationWindow Class Reference

Inheritance diagram for PlotAnimations.AnimationWindow:



Public Member Functions

- `def __init__ (self, parent, datalist, config=None, yvals=None, mode="1D", args, kwargs)`
- `def on_close (self, e)`
- `def update (self, frame_number)`
- `def newplot (self)`
- `def init (self)`
- `def on_play (self, e)`
- `def RefreshPlot (self)`
- `def on_next (self, e)`
- `def on_back (self, e)`
- `def update_framerate (self, e)`

Public Attributes


- `mode`
- `config`
- `datalist`
- `yvals`
- `dim`
- `pos`
- `play`
- `panel`
- `sizer`
- `plot`
- `controlsizer`
- `sb`
- `sbs`
- `frslider`
- `playbutton`
- `nextbutton`
- `backbutton`
- `ctlautoscale`
- `data`
- `xlim`
- `yylim`
- `animation`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/PlotAnimations.py`

3.3 masstools.AutocorrWindow Class Reference

Inheritance diagram for masstools.AutocorrWindow:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, config, massdat)
- def **Go** (self, e)
- def **PickPeaks** (self, e)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)

Public Attributes


- **config**
- **massdat**
- **pnl**
- **vbox**
- **hbox**
- **sb**
- **sbs**
- **plot1**
- **listpanel**
- **corr**
- **corrx**
- **pks2**
- **peaks**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.4 plot1d.BarChart Class Reference

Inheritance diagram for plot1d.BarChart:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **plotrefreshtop** (self, X, Y, titlestr, xstr, ystr, colortab, peakval)
- def **plotadd** (self, X, Y, colval, newlabel)
- def **plotaddlegend** (self)
- def **plotadddot** (self, X, Y, colval, markval)
- def **PaintAgain** (self, bins, kwargs)
- def **PlotClear** (self)
- def **on_save_fig_dialog** (self, evt)
- def **on_save_fig** (self, evt, path, kwargs)

Public Attributes


- **flag**
- **subplot1**
- **labels**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/plot1d.py

3.5 nativez.ColorList Class Reference

Inheritance diagram for nativez.ColorList:



Public Member Functions

- def **__init__** (self, parent)
- def **AddLineEmpty** (self)
- def **AddLine** (self, zoff)
- def **ReturnData** (self)
- def **SuperDelete** (self)

Public Attributes

- **ultimateList**
- **buttontot**
- **colorbox**
- **zoffouts**

3.5.1 Detailed Description

3.5.2 Constructor & Destructor Documentation

3.5.2.1 `def nativez.ColorList.__init__(self, parent)`


Constructor

The documentation for this class was generated from the following file:

- C:/Python/UniDec/nativez.py

3.6 ColorPlot.ColorPlot2D Class Reference

Inheritance diagram for ColorPlot.ColorPlot2D:



Public Member Functions

- `def __init__(self, args, kwargs)`
- `def ColorPlot(self, mzgrid, mzax, dtax, ztab)`
- `def PaintAgain(self, bins, kwargs)`
- `def on_save_fig(self, evt, path, kwargs)`
- `def PlotClear(self)`

Public Attributes


- `flag`
- `data_x_lim`
- `data_y_lim`
- `data_lims`
- `mzgrid`
- `mzlen`
- `dtlen`
- `zlen`
- `ztot`
- `zind`
- `ztab`
- `skew`
- `topcmap`
- `colors`
- `subplot1`
- `norm`
- `cax`
- `cbar`

The documentation for this class was generated from the following file:

- C:/Python/UniDec/ColorPlot.py

3.7 masstools.CorrListCtrl Class Reference

Inheritance diagram for masstools.CorrListCtrl:



Public Member Functions


- def **__init__**
- def **Populate** (self, pks)
- def **Clear** (self)

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.8 masstools.CorrListCtrlPanel Class Reference

Inheritance diagram for masstools.CorrListCtrlPanel:



Public Member Functions

- def **__init__** (self, parent)
- def **OnUseNative** (self, event)

Public Attributes


- **useNative**
- **list**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.9 plot3d.CubePlot Class Reference

Inheritance diagram for plot3d.CubePlot:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **cubeplot** (self, xaxis, yaxis, zaxis, C, C2, C3, xlab, ylab, zlab, config)
- def **on_save_fig_dialog** (self, evt)
- def **on_save_fig** (self, evt, path, kwargs)
- def **PaintAgain** (self, bins)
- def **PlotClear** (self)
- def **MakeIsoMatrices** (self)
- def **IsometricProjection** (self, x, y, z)
- def **Isogrids** (self)
- def **IsoTicks** (self, ax, col)
- def **IsoLabs** (self, xtclab, yticlab, zticlab, ax, col)
- def **IsoLines** (self)

Public Attributes


- **flag**
- **subplot1**
- **alpha**
- **beta**
- **bins**
- **N**
- **cmap**
- **col**
- **xaxis**
- **yaxis**
- **zaxis**
- **xlen**
- **ylen**
- **zlen**
- **xticlab**
- **yticlab**
- **zticlab**
- **xticloc**
- **yticloc**
- **zticloc**
- **cbar**
- **datalims**
- **A**
- **B**
- **P**
- **C**
- **PC**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/plot3d.py

3.10 datacollector.DataCollector Class Reference

Inheritance diagram for datacollector.DataCollector:



Public Member Functions

- def **__init__** (self, parent, title, config=None, pks=None, args, kwargs)
- def **LoadXfromPeaks** (self, e)
- def **on_save** (self, e)
- def **on_load** (self, e)
- def **load** (self, savename)
- def **on_add_x** (self, e)
- def **on_add_y** (self, e)
- def **on_choose_dir** (self, e)
- def **OnMotion** (self, xpos, ypos)
- def **update_get** (self, e)
- def **update_set** (self, e)
- def **data_extract**
- def **on_run**
- def **MakePlot2** (self)
- def **on_save_fig** (self, e)
- def **on_save_figPDF** (self, e)
- def **on_kd_fit** (self, e)
- def **on_animate** (self, e)
- def **on_animate2** (self, e)
- def **on_2dgrid** (self, e)
- def **on_defect** (self, e)
- def **on_autocorr** (self, e)
- def **on_local_path** (self, e)
- def **on_absolute_path** (self, e)
- def **on_export** (self, e)
- def **on_MSMS_norm** (self, e)

Public Attributes

- **directory**
- **config**
- **pks**
- **gridparams**
- **filemenu**
- **menuSave**
- **menuLoad**
- **menuSaveFigPNG**
- **menuSaveFigPDF**
- **toolsmenu**
- **experimentalmenu**
- **menuAnimation**

- menuAnimation2
- menu2dGrid
- menudefect
- menulocalpath
- menuabsolutePath
- menumsmsnorm
- menuautocorr
- menuBar
- panel
- sizer
- inputsizer
- ypanelsizer
- ypanel
- ypanelsizer2
- addybutton
- dirinput
- dirbutton
- xpanel
- xpanelsizer
- addxbutton
- plotwindow
- tab1
- tab2
- tab3
- plot1
- plot2d
- plot4
- plot2
- plotwindow2
- tab12
- tab22
- plot3
- plot3h
- runsizer
- ctldata
- ctlmin
- ctlmax
- ctlnorm
- runbutton
- ctlnorm2
- ctlextract
- ctlwindow
- runsizer2
- ctlprot
- ctllig
- ctlbootstrap
- kdbutton
- ctlprotmodel
- ctlligmodel
- ctlmaxsites
- ctloutliers
- plotsizer
- xvals
- yvals
- range

- **extract**
- **normflag**
- **normflag2**
- **protflag**
- **ligflag**
- **datachoice**
- **numprot**
- **numlig**
- **bootstrap**
- **window**
- **maxsites**
- **extractchoice**
- **savename**
- **localpath**
- **molig**
- **data**
- **var1**
- **grid**
- **xlabel**
- **xcolors**
- **aniwindow**
- **compress**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/datacollector.py

3.11 unidec.DataContainer Class Reference

Public Member Functions

- def [__init__](#) (self)

Public Attributes

- **fitdat**
- **fitdat2d**
- **rawdata**
- **rawdata3**
- **data2**
- **data3**
- **massdat**
- **mzgrid**
- **massgrid**
- **ztab**
- **massccs**
- **ccsz**
- **ccsdata**

3.11.1 Constructor & Destructor Documentation

3.11.1.1 def unidec.DataContainer.__init__(self)


Initialize DataContainer with empty arrays.
:return: None

The documentation for this class was generated from the following file:

- C:/Python/UniDec/unidec.py

3.12 Extract2D.Extract2DPlot Class Reference

Inheritance diagram for Extract2D.Extract2DPlot:



Public Member Functions

- def **__init__** (self, parent, datalist, config=None, yvals=None, directory=None, header=None, params=None, args, kwargs)
- def **modparams** (self)
- def **getfromgui** (self)
- def **makegrid** (self)
- def **extractall** (self)
- def **makeplot** (self)
- def **makeplottotal** (self)
- def **makeplotwap** (self)
- def **on_close** (self, e)
- def **on_back** (self, e)
- def **on_next** (self, e)
- def **on_total** (self, e)
- def **on_wap** (self, e)
- def **on_save_fig** (self, e)
- def **on_save_figPDF** (self, e)

Public Attributes

- **filemenu**
- **menuSaveFigPNG**
- **menuSaveFigPDF**
- **menuBar**
- **config**
- **directory**
- **header**
- **params**
- **datalist**
- **dlen**


- pos
- yvals
- storediscrete
- panel
- sizer
- plot1
- plot2
- controlsizer
- controlsizer1
- ctlm0
- ctlm1
- ctlm2
- ctlm1min
- ctlm1max
- ctlm2min
- ctlm2max
- ctlwindow
- ctlnorm
- controlsizer2
- backbutton
- nextbutton
- totalbutton
- wapbutton
- m0
- m1
- m2
- m1minmax
- m2minmax
- window
- normflag
- grids
- m1range
- m2range
- m2grid
- massgrid
- igrd

The documentation for this class was generated from the following file:

- C:/Python/UniDec/Extract2D.py

3.13 miscwindows.FileNameDialog Class Reference

Inheritance diagram for miscwindows.FileNameDialog:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, config)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)

Public Attributes


- **config**
- **defaultin**
- **defaultout**
- **defaultconf**
- **defaultmassfile**
- **defaulttruncfile**
- **defaultofile**
- **defaultmatchfile**
- **defaultpeaksfile**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **hbox1**
- **inputbox**
- **hbox2**
- **confbox**
- **hbox3**
- **outputbox**
- **hbox4**
- **massbox**
- **hbox5**
- **truncbox**
- **hbox6**
- **obox**
- **hbox7**
- **matchbox**
- **hbox8**
- **peakbox**
- **inboxval**
- **outboxval**
- **confboxval**
- **truncval**
- **massfileval**
- **ofileval**
- **matchfileval**
- **peaksfileval**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/miscwindows.py

3.14 import_wizard.ImportWizard Class Reference

Inheritance diagram for import_wizard.ImportWizard:



Public Member Functions

- def **__init__**
- def **setup_frame** (self)
- def **get_folder_path** (self, evt)
- def **on_folder_path_change** (self, evt)
- def **add_file** (self, evt)
- def **export_then_load** (self, evt)
- def **auto**
- def **export_file** (self, evt)
- def **close** (self, evt)

Public Attributes

- **exedir**
- **Title**
- **file_set**
- **rb**
- **folder_path**
- **my_grid**
- **my_tree**
- **tree**
- **desc**

3.14.1 Member Function Documentation

3.14.1.1 def import_wizard.ImportWizard.export_file (self, evt)

Export import file for later import

3.14.1.2 def import_wizard.ImportWizard.export_then_load (self, evt)

save file then load it

3.14.1.3 def import_wizard.ImportWizard.get_folder_path (self, evt)


Get path to folder and place in txtctrl

The documentation for this class was generated from the following file:

- C:/Python/UniDec/import_wizard.py

3.15 IM_windows.IMToolExtract Class Reference

Inheritance diagram for IM_windows.IMToolExtract:



Public Member Functions

- `def __init__ (self, args, kwargs)`
- `def InitUI (self, massdat, ccscat, mcsgrid, config, pks)`
- `def OnClose (self, e)`
- `def OnCloseCancel (self, e)`
- `def loadpeaks (self, e)`
- `def getfromgui (self, e)`
- `def OnAdd (self, e)`
- `def OnPlot (self, e)`

Public Attributes


- `config`
- `massdat`
- `ccscat`
- `totalgrid`
- `pks`
- `ztab`
- `zstrings`
- `pnl`
- `vbox`
- `sb`
- `sbs`
- `plot1`
- `plot2`
- `plotsizer`
- `sb2`
- `sbs2`
- `gbox1c`
- `ctlzout`
- `masspanel`
- `addbutton`
- `plotbutton`
- `hboxend`
- `okButton`
- `closeButton`
- `zout`
- `zoutgrid`
- `ccsproj`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/IM_windows.py`

3.16 IM_windows.IMTools Class Reference

Inheritance diagram for IM_windows.IMTools:



Public Member Functions

- `def __init__ (self, args, kwargs)`
- `def InitUI (self, data3, config)`
- `def SetupPanel (self)`
- `def OnClose (self, e)`
- `def OnCloseCancel (self, e)`
- `def loadtogui (self, e)`
- `def getfromgui (self, e)`
- `def OnAdd (self, e)`
- `def OnPlot (self, e)`
- `def OnFlipTWave (self, e)`

Public Attributes

- `defaultconfig`
- `config`
- `data3`
- `pnl`
- `pnl2`
- `flag`
- `vbox`
- `sb`
- `sbs`
- `plot`
- `ctlsizer`
- `sb2`
- `sbs2`
- `gbox1c`
- `ctltwave`
- `ctlvolt`
- `ctlpressure`
- `ctltemp`
- `ctlgasmass`
- `ctlto`
- `ctldriftlength`
- `twave`
- `ctltcal1`
- `ctltcal2`
- `ctiledc`
- `vbox2`
- `masspanel`
- `addbutton`

- **plotbutton**
- **hboxend**
- **okButton**
- **closeButton**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/IM_windows.py

3.17 UniFit.KDmodel Class Reference

Public Member Functions

- def **__init__** (self, numtotprot, numtotlig, data, pconc, lconc, nodelist, header, removeoutliers=False, plot1=None, plot2=None, plot3=None, bootnum=1, maxsites=0, kwargs)
- def **FindBestModel** (self, fixedprotmodel, fixedligmodel)
- def **RunKDFit** (self, kwargs)
- def **Return** (self)
- def **SetupModel** (self, kwargs)
- def **ModPaths** (self)
- def **MakeFitGrid** (self)
- def **OutlierTest** (self)
- def **Bootstrap** (self, std, numpts)
- def **RunBootstrap** (self)
- def **GraphPlot**
- def **PlotTrace**
- def **PlotHist**

Public Attributes

- **outlierflag**
- **plot1**
- **plot2**
- **plot3**
- **protflag**
- **ligflag**
- **mode**
- **header**
- **maxsites**
- **kdargs**
- **randfit**
- **bootnum**
- **numtotprot**
- **numtotlig**
- **nprot**
- **nlig**
- **data**
- **fixedligmodel**
- **fixedprotmodel**
- **findmodelflag**
- **plotflag**
- **ligmodel**
- **fit**

- **degenkds**
- **residuals**
- **stdres**
- **error**
- **reactions**
- **nodes**
- **kdmap**
- **numkd**
- **graph**
- **ugraph**
- **nodenames**
- **fitgrid**
- **stddevs**
- **meanfit**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/UniFit.py

3.18 UniFit.kdstruct Class Reference

Public Member Functions

- **def __init__** (self)

Public Attributes


- **pconc**
- **lconc**
- **ureact**
- **nprottab**
- **nligtab**
- **paths**
- **pfrees**
- **lfrees**
- **weights**
- **odelist**
- **kds**
- **degen**
- **maxsites**
- **nfactors**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/UniFit.py

3.19 datacollector.ListCtrlPanel Class Reference

Inheritance diagram for datacollector.ListCtrlPanel:



Public Member Functions

- `def __init__ (self, parent, type="X", size=(200, 400))`
- `def OnUseNative (self, event)`
- `def OnRightClick (self, event)`
- `def OnPopupOne (self, event)`
- `def OnPopupTwo (self, event)`
- `def OnPopupThree (self, event)`
- `def OnPopupFour (self, event)`

Public Attributes


- `useNative`
- `list`
- `popupID1`
- `popupID2`
- `popupID3`
- `popupID4`
- `selection`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/datacollector.py`

3.20 mainwindow.Mainwindow Class Reference

Inheritance diagram for mainwindow.Mainwindow:



Public Member Functions

- `def __init__ (self, parent, title, config)`
- `def setup_shortcuts (self)`
- `def setup_menu (self)`
- `def on_defaults (self, e)`
- `def setup_main_panel (self)`

- def [setup_tool_tips](#) (self)
- def [import_config_to_gui](#) (self)
- def [export_gui_to_config](#) (self)
- def [menu_401_403](#) (self, event)
- def [clear_all_plots](#)
- def [on_motion](#) (self, xpos, ypos)
- def [on_flip_mode](#) (self, e)
- def [on_flip_tabbed](#) (self, e)
- def [on_flip_twave](#) (self, e)
- def [on_about](#) (self, e)
- def [on_exit](#) (self, e)
- def [save_all_figures](#) (self, extension, extension2="", e=0, header=None, kwargs)
- def [on_save_figure_eps](#) (self, e)
- def [on_save_figure_png](#) (self, e, kwargs)
- def [on_save_figure_pdf](#) (self, e)
- def [on_save_figure_dialog](#) (self, e)
- def [shrink_figure](#) (self, plot)
- def [shrink_all_figures](#) (self)
- def [on_save_figure_small](#) (self, e)
- def [on_check_manual](#) (self, e)
- def [on_mass_list](#) (self, e)

Public Attributes

- **pres**
- **config**
- **open_bmp**
- **next_bmp**
- **report_bmp**
- **A_bmp**
- **ud_bmp**
- **system**
- **displaysize**
- **tabbed**
- **imflag**
- **twave**
- **filemenu**
- **toolsmenu**
- **analysismenu**
- **advancedmenu**
- **experimentalmenu**
- **menuOpen**
- **menuOpenRaw**
- **menuLoadState**
- **menuSaveState**
- **menuLoad**
- **menuLoadDefault**
- **menuSaveDefault**
- **defaultmenu**
- **menuDefault1**
- **menuDefault2**
- **menufigdialog**
- **figmenu**
- **menuSaveFigure0**

- menuSaveFigure1s
- menuSaveFigure1
- menuSaveFigure2
- menuSaveFigure4
- menuAbout
- menuExit
- menuBatch
- menuBatch2
- menuBatchRaw
- menuImportWizard
- menuWidth
- menuManualFile
- menuMassFile
- menuPlotZ
- menucollect
- menuExport
- menuFitNorm
- menukendrick
- menu2Dgrid
- menuintegrate
- menumatch
- menucom
- menucolor1d
- menuoffset
- menuimtools
- menuimtools2
- menunativeccs
- menuReset
- menuUnidecPath
- menuFileName
- menuflipmode
- menufliptabbed
- Tweet
- menuAdditionalParameters
- menuDeisotope
- menuCrossValidate
- menusuperbatch
- menumassprocess
- menupastespectrum
- menuerrors
- menuBar
- plot1
- plot2
- plot3
- plot4
- plot5
- plot6
- plot1im
- plot1fit
- plot2ccs
- plot5mccs
- plot5ccsz
- plot3color
- plot9
- plot10

- peakpanel
- openbutton
- procbutton
- udbutton
- ppbutton
- autobutton
- ctlminmz
- ctlmaxmz
- ctlmindt
- ctlmaxdt
- ctlsmoothdt
- ctlsubbuffdt
- subtypectl
- dataprepbutton
- ctlbuff
- ctlsmooth
- ctlbinsize
- ctlpusher
- ctlinththresh
- ctladductmass
- ctllacelvolt
- ctlbintype
- ctlconvertflag
- ciltwave
- ctlvolt
- ctlpressure
- ctltemp
- ctlgasmass
- ctlto
- ctldriflength
- ciltcal1
- ciltcal2
- ctledc
- ctlistartz
- ctlendz
- ctlmasslb
- ctlmassub
- ctlccslb
- ctlcsub
- ctlcsubins
- ctldtsig
- ctlmassbins
- ctlmzsig
- ctlipsfun
- rununidec
- ctlzsig
- ctlmolig
- ctlmsig
- ctlcsg
- ctlnumit
- ctlpoolflag
- ctlisotopemode
- ctlmanualassign
- ctlmasslistflag
- ctltabsig

- `ctlminnativez`
- `ctlmaxnativez`
- `ctlnativeccslb`
- `ctlnativeccsub`
- `ctlwindow`
- `ctlthresh`
- `ctlnorm`
- `plotbutton`
- `plotbutton2`
- `ctl2dcm`
- `ctlpeakcm`
- `ctldiscrete`
- `ctlpublicationmode`
- `ctlrawflag`
- `ctlthresh2`
- `ctlsep`
- `ctlintlb`
- `ctlintub`
- `replotbutton`
- `compositebutton`
- `cubeplotbutton`
- `pngs`
- `figsize`
- `rect`

3.20.1 Detailed Description

Main UniDec GUI Window.

3.20.2 Constructor & Destructor Documentation

3.20.2.1 `def mainwindow.Mainwindow.__init__(self, parent, title, config)`

initialize window and feed in links to presenter and config.

```
:param parent: GUniDec Presenter -> self.pres
:param title: Window title (string)
:param config: UniDecConfig object ->self.config
:return: None
```

3.20.3 Member Function Documentation

3.20.3.1 `def mainwindow.Mainwindow.clear_all_plots(self, flag = 0)`

Clear All Plots
:return: None

3.20.3.2 `def mainwindow.Mainwindow.export_gui_to_config(self)`

Exports parameters from the GUI to the config object.
:return: None

3.20.3.3 def mainwindow.Mainwindow.import_config_to_gui(self)

Imports parameters from the config object to the GUI.
:return: None

3.20.3.4 def mainwindow.Mainwindow.menu_401_403(self, event)

Menu function to adjust the UniDec core function (agreesiveflag).
:param event: wx Event
:return: None

3.20.3.5 def mainwindow.Mainwindow.on_about(self, e)

Displays message about program
:param e:
:return:

3.20.3.6 def mainwindow.Mainwindow.on_check_manual(self, e)

Checks the configuration to see if values for manual mode are set. If they are not, it opens the window to set the manual assignments.
:param e: Dummy wx event passed on.
:return: None

3.20.3.7 def mainwindow.Mainwindow.on_defaults(self, e)

Resets the configuration to a default predefined in the unidecstructure file.
:param e: Menu event
:return: None

3.20.3.8 def mainwindow.Mainwindow.on_exit(self, e)

Exit the Program
:param e: Dummy wx event
:return: None

3.20.3.9 def mainwindow.Mainwindow.on_flip_mode(self, e)

Flips between MS and IM-MS mode
:param e: wx event or anything (will flip if not 0)
:return: None

3.20.3.10 def mainwindow.Mainwindow.on_flip_tabbed(self, e)

Flips between tabbed plots and a single window of plots
:param e: wx Event or anything (will flip if not 0)
:return: None

3.20.3.11 def mainwindow.Mainwindow.on_flip_twave(self, e)

Flips between T-Wave and Linear IM-MS
:param e: wx Event or anything (will get value from Selection if not 0)
:return: None

3.20.3.12 def mainwindow.Mainwindow.on_mass_list (self, e)

Checks the configuration to see if values for the mass list are set. If they are not, it opens the window to set the mass list.

```
:param e: Dummy wx event passed on.
:return: None
```

3.20.3.13 def mainwindow.Mainwindow.on_motion (self, xpos, ypos)

Triggered by pubsub from plot windows. Reports text in Status Bar.

```
:param xpos: x position fed from event
:param ypos: y position fed from event
:return: None
```

3.20.3.14 def mainwindow.Mainwindow.on_save_figure_dialog (self, e)

Open dialog box to set the parameters for figure type, size, and path to save.

```
:param e: Dummy wx event
:return: None
```

3.20.3.15 def mainwindow.Mainwindow.on_save_figure_eps (self, e)

Save all figures as EPS

```
:param e: Dummy wx event
:return: None
```

3.20.3.16 def mainwindow.Mainwindow.on_save_figure_pdf (self, e)

Saves all figures as PDF

```
:param e: Dummy wx event
:return: None
```

3.20.3.17 def mainwindow.Mainwindow.on_save_figure_png (self, e, kwargs)

Save all figures as PNG

```
:param e: Dummy wx event
:param kwargs: keywords to pass to matplotlib savefig
:return: None
```

3.20.3.18 def mainwindow.Mainwindow.on_save_figure_small (self, e)

Preset to shrink figures to 4.5 in by 3 in and save as PDF.

```
:param e: Dummy wx event
:return: None
```

3.20.3.19 def mainwindow.Mainwindow.save_all_figures (self, extension, extension2 = ' ', e = 0, header = None, kwargs)

Save All of the Figures. Will name as header+extension2+_FigureX.+extension

```
:param extension: Figure type (pdf, eps, png). Anything accepted by matplotlib
:param extension2: Additional text to include in the figure header.
:param e: Dummy wx Event
:param header: Option to add different header. Default of none yields self.outfname as the path header
:param kwargs: Any keywords to pass to the matplotlib savefig command such as Transparent or DPI
:return: figureflags, files (the figures that were successfully saved and the files that they were saved to)
```

3.20.3.20 `def mainwindow.Mainwindow.setup_main_panel (self)`

Lays Out Main Panel. Binds some functions to presenter.
:return: None

3.20.3.21 `def mainwindow.Mainwindow.setup_menu (self)`

Sets menu and binds menu objects to functions in presenter and window
:return: None

3.20.3.22 `def mainwindow.Mainwindow.setup_shortcuts (self)`

Setup shortcuts in GUI. Binds key combinations to functions in presenter (self.pres)
:return: None

3.20.3.23 `def mainwindow.Mainwindow.setup_tool_tips (self)`

Sets Tool Tips for items on the Main Panel
:return: None

3.20.3.24 `def mainwindow.Mainwindow.shrink_all_figures (self)`

Shrinks all figures to the size specified in self. figsize
:return: A list of plot objects that we shrunk

3.20.3.25 `def mainwindow.Mainwindow.shrink_figure (self, plot)`


Automatically shrinks the plot to a figure size in inches set in self. figsize.
:param plot: Plot object to shrink
:return: None

The documentation for this class was generated from the following file:

- C:/Python/UniDec/mainwindow.py

3.21 masstools.ManualSelection Class Reference

Inheritance diagram for masstools.ManualSelection:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, config, data)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)
- def **OnClear** (self, e)
- def **OnAdd** (self, e)
- def **OnAddFromPlot** (self, e)
- def **OnPlot** (self, e)
- def **OnImport** (self, e)

Public Attributes

- **data**
- **config**
- **defaulttrunc**
- **newtrunc**
- **pnl**
- **vbox**
- **vbox2**
- **plot1**
- **hbox**
- **sb**
- **sbs**
- **importbutton**
- **clearbutt**
- **addbutton**
- **addbutton2**
- **plotbutton**
- **sb2**
- **sbs2**
- **masslistbox**
- **truncfilename**
- **importtrunc**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.22 MassModelFitter.mass Class Reference

Public Member Functions

- def **__init__** (self, args)
- def **LoadGuesses** (self, mass0)
- def **MakePeakList** (self)
- def **f0** (self, array, peaks, run)
- def **f0s** (self, array, peaks, run)
- def **ComparePeaks** (self, peaks)

Public Attributes


- **mass0**
- **mass1**
- **mass2**
- **dist1**
- **dist1limits**
- **dist2**
- **dist2limits**
- **adduction**
- **psfun**
- **resolution**
- **thresh**
- **intensity**
- **num1**
- **num2**
- **int1**
- **int2**
- **num2grid**
- **massgrid**
- **int2grid**
- **intgrid**
- **masspeaks**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/MassModelFitter.py

3.23 MassDefects.MassDefectWindow Class Reference

Inheritance diagram for MassDefects.MassDefectWindow:



Public Member Functions

- def **__init__** (self, parent, datalist, config=None, yvals=None, pks=None, value=None, dir=None, show=True, args, kwargs)
- def **OnClose** (self, e)
- def **getfromgui** (self)
- def **makegrid** (self)
- def **extractall** (self)
- def **makeplot** (self)
- def **makeplottotal** (self)
- def **on_back** (self, e)
- def **on_next** (self, e)
- def **on_total** (self, e)
- def **on_peaks** (self, e)
- def **on_save_fig** (self, e)
- def **on_save_figPDF** (self, e)
- def **on_add_line** (self, e)

Public Attributes

- **directory**
- **parent**
- **filemenu**
- **menuSaveFigPNG**
- **menuSaveFigPDF**
- **plotmenu**
- **menuaddline**
- **menuBar**
- **config**
- **datalist**
- **pos**
- **yvals**
- **panel**
- **sizer**
- **plot1**
- **plot2**
- **controlsizer**
- **ctlm0**
- **ctlwindow**
- **controlsizer2**
- **backbutton**
- **nextbutton**
- **totalbutton**
- **pks**
- **peaksbutton**
- **radiobox**
- **radiobox2**
- **radiobox3**
- **ylab**
- **m0**
- **window**
- **ktype**
- **factor**
- **xlab**
- **xaxis**
- **kmass**
- **nominalkmass**
- **kmdefectexact**
- **defects**
- **nominal**
- **m2grid**
- **igrid**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/MassDefects.py

3.24 MassFitter.MassFitter Class Reference

Public Member Functions

- **def __init__** (self, massdat, finarray, psfun, args)
- **def Fit** (self, args)

Public Attributes


- **massdat**
- **finarray**
- **psfun**
- **initguess**
- **fit**
- **fitdat**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/MassFitter.py

3.25 masstools.MassSelection Class Reference

Inheritance diagram for masstools.MassSelection:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI**
- def **OnSim** (self, e)
- def **OnPlot** (self, e)
- def **OnRightClick** (self, event)
- def **OnPopupOne** (self, event)
- def **OnRightClick2** (self, event)
- def **OnPopup2** (self, event)
- def **OnMatchI** (self, e)
- def **OnMatchAll** (self, e)
- def **Match** (self)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)
- def **OnPopulateButton** (self, e)
- def **OnPopulateButton2** (self, e)
- def **OnPopulateButton3** (self, e)
- def **OnClear** (self, e)
- def **OnAdd** (self, e)
- def **OnAdd2** (self, e)
- def **OnClear2** (self, e)
- def **OnImport** (self, e)
- def **OnImport2** (self, e)

Public Attributes

- **massdat**
- **config**
- **defaultmasslist**
- **defaultoligolist**
- **defaultmatchlist**
- **newmasslist**
- **oligos**
- **matchlist**
- **newmatchlist**
- **pks**
- **pnl**
- **vbox**
- **hbox**
- **sb**
- **sbs**
- **peakpop**
- **importbutton**
- **oligopop**
- **oligopop2**
- **clearbutt**
- **addbutton**
- **simbutton**
- **masslistbox**
- **sb2**
- **sbs2**
- **clearbutt2**
- **addbutton2**
- **importbutton2**
- **plotbutton**
- **buttonbox**
- **hbox3**
- **oligomerlistbox**
- **sb4**
- **sbs4**
- **matchlbutt**
- **matchAllbutt**
- **matchlistbox**
- **hbox2**
- **newpeaks**
- **popupID1**
- **selection**
- **popupID2**
- **oligomasslist**
- **oligonames**
- **matches**
- **error**
- **peaks**
- **names**
- **oligoshort**
- **mfilename**
- **importmass**
- **ofilename**
- **importolig**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.26 mzMLImporter.mzMLImporter Class Reference

Public Member Functions

- def **__init__** (self, path, args, kwargs)
- def **GetData** (self)

Public Attributes


- **msrun**
- **data**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/mzMLImporter.py

3.27 nativez.NativeZ Class Reference

Inheritance diagram for nativez.NativeZ:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, xvals, yvals, zdat, config, pks)
- def **OnReplot** (self, e)
- def **PlotZoffs** (self)
- def **MakePlot7** (self, e)
- def **on_edit** (self, event)
- def **on_delete** (self, event)
- def **fit** (self, e)
- def **update** (self, e)
- def **OnReset** (self, e)
- def **UpdateList** (self)
- def **PopulateList** (self, e)
- def **onadd** (self, e)
- def **MakeFArray** (self, min, max)
- def **fastextract** (self, f, width)
- def **Extract** (self, e)
- def **GetMaxima** (self)
- def **PeakExtract** (self)
- def **SaveFig** (self, e)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)

Public Attributes


- **config**
- **pks**
- **xlen**
- **ylen**
- **xvals**
- **yvals**
- **newgrid**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **hbox0**
- **plot1**
- **plot2**
- **plot3**
- **plot4**
- **plot5**
- **plot6**
- **plot7**
- **hbox1**
- **addbutton**
- **fitbutton**
- **resetbutton**
- **extractbutton**
- **massoffset**
- **ctlfilt**
- **savefigbutt**
- **replotbutton**
- **hbox2**
- **list**
- **hbox3**
- **zoffs**
- **eshape**
- **offsetgrid**
- **ftot**
- **f**
- **offsetvalue**
- **extracts**
- **maxes**
- **peakextracts**
- **peakextractsarea**
- **zoffouts**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/nativez.py

3.28 datacollector.NetworkFrame Class Reference

Inheritance diagram for datacollector.NetworkFrame:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **on_save_fig** (self, evt, path, kwargs)
- def **Clear** (self)

Public Attributes

- **axes**
- **flag**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/datacollector.py

3.29 peakstructure.Peak Class Reference

Public Member Functions

- def **__init__** (self)

Public Attributes

- **mass**
- **height**
- **ccs**
- **area**
- **color**
- **label**
- **marker**
- **textmarker**
- **ignore**
- **match**
- **matcherror**
- **integral**
- **integralrange**
- **mztab**
- **mztab2**
- **stickdat**

- **kendricknum**
- **kendrickdefect**
- **kmass**
- **score**
- **corrint**
- **correrr**
- **mztabi**
- **massavg**
- **masserr**
- **tval**
- **peakmasses**
- **fitmassavg**
- **fitmasserr**
- **fitarea**
- **fitareaerr**

3.29.1 Detailed Description

Class for a single peak. Contains all key parameters for describing and plotting the peak.

3.29.2 Constructor & Destructor Documentation

3.29.2.1 `def peakstructure.Peak.__init__(self)`

Initialize all parameters for the peak to defaults

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peakstructure.py

3.30 peakstructure.Peaks Class Reference

Public Member Functions

- `def __init__(self)`
- `def add_peaks(self, parray)`
- `def default_params`
- `def get_mass_defects`
- `def score_peaks`

Public Attributes

- **peaks**
- **plen**
- **changed**
- **masses**
- **convolved**
- **composite**
- **peakcolors**
- **markers**
- **colormap**
- **textmarkers**
- **marklen**

3.30.1 Detailed Description

Class containing all useful data about peaks.

The peaks themselves are of the Peak class and contained within the `self.peaks` list.

3.30.2 Constructor & Destructor Documentation

3.30.2.1 `def peakstructure.Peaks.__init__(self)`

Initialize Peaks class and set empty values
:return: None

3.30.3 Member Function Documentation

3.30.3.1 `def peakstructure.Peaks.add_peaks(self, parray)`

Create peak objects from an array
:param parray: N x 2 array containing (mass, height) of each peak.
:return: None

3.30.3.2 `def peakstructure.Peaks.default_params(self, cmap = "rainbow")`

Set default parameters for peaks, such as color, label, and marker
:param cmap: Colormap from matplotlib.cm
:return: None

3.30.3.3 `def peakstructure.Peaks.get_mass_defects(self, kendrickmass, mode = 0)`

Get the mass defects and mass number for each peak
:param kendrickmass: Kendrick reference mass
:param mode: Select range of defects 0=(0,1), 1=(-0.5,0.5)
:return: None

3.30.3.4 `def peakstructure.Peaks.score_peaks(self, thresh = 0, ci = 0.99)`


For each peak, assign a score of the fractional number of charge states observed.
:param thresh: Optional threshold to define a noise floor.
:return: None

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peakstructure.py

3.31 peakwidthtools.PeakTools1d Class Reference

Inheritance diagram for peakwidthtools.PeakTools1d:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, config, data)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)
- def **OnReset** (self, e)
- def **OnCenter** (self, e)
- def **OnPlot** (self, e)
- def **OnFit** (self, e)

Public Attributes


- **config**
- **data**
- **topmax**
- **topmin**
- **mzsig**
- **length**
- **psfun**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **plot1**
- **hbox11**
- **hbox10**
- **centerbutton**
- **hbox9**
- **ctlpsfun**
- **fitbutton**
- **plotbutton**
- **hbox8**
- **ctlmzsig**
- **hbox6**
- **errorbox**
- **hbox7**
- **resbox**
- **centdat**
- **mzsig2**
- **mid**
- **scale**
- **psguess**
- **res**
- **cent**
- **cent2**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peakwidthtools.py

3.32 peakwidthtools.PeakTools2d Class Reference

Inheritance diagram for peakwidthtools.PeakTools2d:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **InitUI** (self, data3, data2, config)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)
- def **OnFlip** (self, e)
- def **OnReset** (self, e)
- def **OnCenter** (self, e)
- def **OnPlot** (self, e)
- def **OnFit** (self, e)

Public Attributes

- **mz**
- **dt**
- **int**
- **C**
- **data**
- **topdata**
- **topmax**
- **topmin**
- **config**
- **mzsig**
- **dtsig**
- **length**
- **psfun**
- **fitmid**
- **fitsig**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **plot1**
- **hbox11**
- **hbox10**
- **centerbutton**
- **flipbutton**
- **hboxsig**
- **outmzsig**
- **hboxsig2**
- **outdtsig**
- **hbox9**


- **ctlpfun**
- **fitbutton**
- **plotbutton**
- **hbox8**
- **ctlmzsig**
- **hbox6**
- **errorbox**
- **flag**
- **pos**
- **intdt**
- **data2**
- **centdat**
- **mzsig2**
- **mid**
- **scale**
- **psguess**
- **cent**
- **cent2**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peakwidthtools.py

3.33 twitter_interface.PinWindow Class Reference

Inheritance diagram for twitter_interface.PinWindow:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)

Public Attributes


- **pnl**
- **vbox**
- **sbs**
- **inputbox**
- **pin**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/twitter_interface.py

3.34 plot1d.Plot1d Class Reference

Inheritance diagram for plot1d.Plot1d:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **plotrefreshtop** (self, X, Y, titlestr, xstr, ystr, lab, config, args, kwargs)
- def **plotrefreshtopspan** (self, X, Y, titlestr, xstr, ystr, kwargs)
- def **plotrefreshtopbox** (self, X, Y, titlestr, xstr, ystr, args, kwargs)
- def **plotadd** (self, X, Y, colval, newlabel)
- def **plotaddspan** (self, X, Y, colval, newlabel)
- def **plotadddot** (self, X, Y, colval, markval)
- def **plotaddlegend**
- def **addtext** (self, txt, x, y, kwargs)
- def **filledplot** (self, x, y, color)
- def **textremove** (self)
- def **PaintAgain** (self, bins, kwargs)
- def **on_save_fig_dialog** (self, evt)
- def **on_save_fig** (self, evt, path, kwargs)
- def **Histogram**
- def **PlotClear** (self, args)

Public Attributes


- **flag**
- **text**
- **lines**
- **kda**
- **subplot1**
- **X**
- **labels**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/plot1d.py

3.35 plot2d.Plot2d Class Reference

Inheritance diagram for plot2d.Plot2d:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **contourplot**
- def **on_save_fig_dialog** (self, evt)
- def **on_save_fig** (self, evt, path, kwargs)
- def **PaintAgain** (self, bins, kwargs)
- def **PlotClear** (self)

Public Attributes


- **flag**
- **subplot1**
- **ticcol**
- **norm**
- **cbar**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/plot2d.py

3.36 plot2d.Plot2dMass Class Reference

Inheritance diagram for plot2d.Plot2dMass:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **contourplot**
- def **on_save_fig_dialog** (self, evt)
- def **on_save_fig** (self, evt, path, kwargs)
- def **PaintAgain** (self, bins, kwargs)
- def **PlotClear** (self)
- def **PlotNativeZ**

Public Attributes


- **flag**
- **nativez**
- **kda**
- **xvals**
- **yvals**
- **subplot1**
- **ticcol**
- **norm**
- **cbar**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/plot2d.py

3.37 PlottingWindow.PlottingWindow Class Reference

Inheritance diagram for PlottingWindow.PlottingWindow:



Public Member Functions

- **def __init__** (self, args, kwargs)
- **def set_link** (self, link)
- **def contour** (self)
- **def save_contour** (self, path, kwargs)
- **def repaint** (self)
- **def set_color**
- **def size_handler** (self, args, kwargs)
- **def setup_zoom**

Public Attributes


- **figure**
- **int**
- **smash**
- **link**
- **canvas**
- **resize**
- **cbar**
- **colors**
- **index2color**
- **cm_dict**
- **zoom**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/PlottingWindow.py

3.38 miscwindows.SaveFigureDialog Class Reference

Inheritance diagram for miscwindows.SaveFigureDialog:



Public Member Functions

- `def __init__ (self, args, kwargs)`
- `def InitUI (self, config)`
- `def OnClose (self, e)`
- `def OnCloseCancel (self, e)`
- `def on_choose_dir (self, e)`

Public Attributes


- `config`
- `directory`
- `header`
- `extension`
- `transparent`
- `rect`
- `figsize`
- `pnl`
- `vbox`
- `sb`
- `sbs`
- `hbox`
- `dirinput`
- `dirbutton`
- `hbox5`
- `headerbox`
- `hbox6`
- `extbox`
- `tbox`
- `hbox7`
- `widebox`
- `tallbox`
- `hbox8`
- `b1`
- `b2`
- `b3`
- `b4`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/miscwindows.py`

3.39 GUniDec.Shell Class Reference

Inheritance diagram for GUniDec.Shell:



Public Member Functions

- `def __init__ (self, args, kwargs)`

Public Attributes


- **shell**
- **shellwindow**

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/GUniDec.py`

3.40 miscwindows.SingleInputDialog Class Reference

Inheritance diagram for miscwindows.SingleInputDialog:



Public Member Functions

- `def __init__ (self, args, kwargs)`
- `def InitUI`
- `def OnClose (self, e)`

Public Attributes


- **pnl**
- **vbox**
- **hbox**
- **inputbox**
- **value**

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/miscwindows.py`

3.41 peaklistsort.TestListCtrl Class Reference

Inheritance diagram for peaklistsort.TestListCtrl:



Public Member Functions


- def **__init__**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peaklistsort.py

3.42 masstools.TestListCtrl Class Reference

Inheritance diagram for masstools.TestListCtrl:



Public Member Functions

- def **__init__**
- def **Populate** (self, listctrlldata)
- def **Clear** (self)
- def **AddLine** (self)
- def **GetList** (self)

Public Attributes


- **currentItem**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.43 masstools.TestListCtrl2 Class Reference

Inheritance diagram for masstools.TestListCtrl2:



Public Member Functions

- def **__init__**
- def **Clear** (self)
- def **AddLine** (self)
- def **Populate** (self, data)
- def **GetList** (self)

Public Attributes


- **index**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.44 masstools.TestListCtrl3 Class Reference

Inheritance diagram for masstools.TestListCtrl3:



Public Member Functions

- def **__init__**
- def **Clear** (self)
- def **AddLine** (self, line=[0])
- def **Populate**
- def **GetList** (self)

Public Attributes


- **imflag**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.45 IM_windows.TestListCtrl4 Class Reference

Inheritance diagram for IM_windows.TestListCtrl4:



Public Member Functions


- def **__init__**
- def **Clear** (self)
- def **AddLine**
- def **Populate**
- def **GetList** (self)

The documentation for this class was generated from the following file:

- C:/Python/UniDec/IM_windows.py

3.46 masstools.TestListCtrlMatch Class Reference

Inheritance diagram for masstools.TestListCtrlMatch:



Public Member Functions

- def **__init__**
- def **Clear** (self)
- def **Populate** (self, data1, data2, data3, data4)
- def **GetList** (self)

Public Attributes


- **index**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.47 masstools.TestListCtrlPanel Class Reference

Inheritance diagram for masstools.TestListCtrlPanel:



Public Member Functions

- def **__init__** (self, parent)
- def **OnUseNative** (self, event)

Public Attributes


- **useNative**
- **list**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.48 peaklistsort.TestListCtrlPanel Class Reference

Inheritance diagram for peaklistsort.TestListCtrlPanel:



Public Member Functions

- def **__init__** (self, parent)
- def **Clear** (self)
- def **AddData**
- def **GetListCtrl** (self)
- def **OnColClick** (self, event)
- def **OnRightClick** (self, event)
- def **OnItemSelected** (self, event)
- def **OnPopupOne** (self, event)
- def **OnPopupTwo** (self, event)
- def **OnPopupThree** (self, event)
- def **OnPopupFour** (self, event)
- def **OnPopupFive** (self, event)

Public Attributes


- **index**
- **list_ctrl**
- **EVT_DELETE_SELECTION_2**
- **EVT_CHARGE_STATE**
- **remove**
- **pks**
- **popupID1**
- **popupID2**
- **popupID3**
- **popupID4**
- **popupID5**
- **currentItem**
- **selection**
- **selection2**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/peaklistsort.py

3.49 masstools.TestListCtrlPanel2 Class Reference

Inheritance diagram for masstools.TestListCtrlPanel2:



Public Member Functions

- def **__init__** (self, parent)
- def **OnUseNative** (self, event)

Public Attributes


- **useNative**
- **list**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/masstools.py

3.50 masstools.TestListCtrlPanel3 Class Reference

Inheritance diagram for masstools.TestListCtrlPanel3:



Public Member Functions

- `def __init__`
- `def OnRightClick (self, e)`
- `def OnPopupOne (self, e)`
- `def OnUseNative (self, event)`

Public Attributes


- `useNative`
- `list`
- `popupID1`
- `selection`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/masstools.py`

3.51 IM_windows.TestListCtrlPanel4 Class Reference

Inheritance diagram for IM_windows.TestListCtrlPanel4:



Public Member Functions

- `def __init__ (self, parent, size=(200, 200))`
- `def OnUseNative (self, event)`

Public Attributes


- `useNative`
- `list`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/IM_windows.py`

3.52 masstools.TestListCtrlPanelMatch Class Reference

Inheritance diagram for masstools.TestListCtrlPanelMatch:



Public Member Functions

- `def __init__ (self, parent)`
- `def OnUseNative (self, event)`

Public Attributes


- `useNative`
- `list`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/masstools.py`

3.53 import_wizard_treectrl.TreeCtrl Class Reference

Inheritance diagram for import_wizard_treectrl.TreeCtrl:



Public Member Functions


- `def __init__ (self, parent, id, pos, size, style)`

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/import_wizard_treectrl.py`

3.54 import_wizard_treectrl.TreeCtrlPanel Class Reference

Inheritance diagram for import_wizard_treectrl.TreeCtrlPanel:



Public Member Functions

- `def __init__ (self, parent, link)`
- `def populate_tree (self)`
- `def add_root`
- `def add_children`
- `def raw_file_info (self, path)`
- `def OnSize (self, event)`
- `def on_selected_changed (self, event)`
- `def on_activate (self, event)`
- `def on_item_expanded (self, event)`

Public Attributes

- `link`
- `tree`
- `fldridx`
- `fldropenidx`
- `fileidx`
- `il`
- `path`
- `raw`
- `item`

3.54.1 Member Function Documentation

3.54.1.1 `def import_wizard_treectrl.TreeCtrlPanel.add_children (self, parent, path, depth_limit = 1)`

Recursively adds children up to the depth specified

3.54.1.2 `def import_wizard_treectrl.TreeCtrlPanel.add_root (self, name, data = None)`


Add root folder to tree
returns root object

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/import_wizard_treectrl.py`

3.55 twitter_interface.TwitterWindow Class Reference

Inheritance diagram for `twitter_interface.TwitterWindow`:



Public Member Functions

- def **__init__** (self, args, kwargs)
- def **LoadScreenName** (self)
- def **OnClose** (self, e)
- def **OnCloseCancel** (self, e)
- def **OnPreview** (self, e)
- def **OnCharacterCount** (self, e)
- def **OnLaunchWeb** (self, e)
- def **Tweet** (self, e)

Public Attributes


- **pngs**
- **codes**
- **imflag**
- **previewsize**
- **APP_KEY**
- **APP_SECRET**
- **pnl**
- **vbox**
- **sb**
- **sbs**
- **hbox1**
- **loginbutton**
- **userbox**
- **hbox2**
- **inputbox2**
- **countbox**
- **hbox3**
- **imagechoice**
- **previewbutton**
- **hbox4**
- **emptyimg**
- **imageCtrl**
- **hbox5**
- **tweetbutton**
- **OAUTH_TOKEN**
- **OAUTH_TOKEN_SECRET**
- **screen_name**
- **imageFile**
- **url**
- **pinwindow**
- **pin**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/twitter_interface.py

3.56 unidec.UniDec Class Reference

Inheritance diagram for unidec.UniDec:



Public Member Functions

- def `__init__` (self)
- def `initialize` (self)
- def `reset_config` (self)
- def `load_config` (self, f_name)
- def `export_config`
- def `open_file` (self, file_name, file_directory, kwargs)
- def `raw_process`
- def `process_data` (self, kwargs)
- def `run_unidec`
- def `unidec_imports`
- def `pick_peaks` (self)
- def `convolve_peaks` (self)
- def `check_badness` (self)
- def `autocorrelation`
- def `kendrick_peaks`
- def `kendrick_continuous`
- def `save_default` (self)
- def `load_default` (self)
- def `mass_grid_to_f_grid` (self)
- def `integrate`
- def `autointegrate`
- def `export_params` (self, e)
- def `process_mass_data` (self)
- def `center_of_mass`
- def `fit_all_masses` (self)
- def `get_charge_peaks` (self)
- def `load_state` (self, load_path)
- def `cross_validate`
- def `normalize_peaks` (self)
- def `align_peaks`
- def `correlate_intensities` (self, pmasses=None, x_range=None, window=None, ci=0.99, kwargs)
- def `get_peaks_scores` (self, window=None, x_range=None, ci=0.99, kwargs)
- def `fit_isolated_peaks` (self, pmasses=None, x_range=None, window=None, norm=False, plot_fits=False, kwargs)
- def `get_errors` (self, kwargs)
- def `open_test_spectrum` (self, masslist=None, n=1, kwargs)
- def `make_plot`

Public Attributes

- **config**
- **data**
- **pks**
- **autopeaks**
- **peakparams**
- **massfit**
- **massfitdat**
- **errorgrid**

3.56.1 Constructor & Destructor Documentation

3.56.1.1 `def unidec.UniDec.__init__(self)`

UniDec Engine

Consists of three main subclasses: Config, DataContainer, Peaks

:return: None

3.56.2 Member Function Documentation

3.56.2.1 `def unidec.UniDec.autocorrelation(self, massdat=None)`

Performs autocorrelation on mass data. Result is stored as `self.data.autocorr`.

Picks peaks greater than 0 using peak detection parameters in config file.

Peaks are stored as a peak structure at `self.autopeaks`

:param massdat: Data on which to run autocorrelation. Default is None, in which case `self.data.massdat` is used

:return: float. First peak in autocorrelation.

3.56.2.2 `def unidec.UniDec.autointegrate(self, ztab=None)`

Perform automatic integration of peaks.

If `self.config.integrateup` is empty, the upperbound becomes `self.config.peakwindow`.

If `self.config.integratelb` is empty, the lowerbound becomes `-self.config.peakwindow`.

Integral range for each peak is set to `peak.integralrange`.

Integral value is set to `peak.integral`.

If `ztab` parameter is set to a list of charge states, it will integrate the mass vs charge grid at each individual charge state. Otherwise, this is ignored.

:param ztab: List of charge states (default = None)

:return: zarea: P x Z array where P is the number of peaks and Z is the number of charge states.

Each value of the array is the integral of peak P at charge state Z.

3.56.2.3 `def unidec.UniDec.center_of_mass(self, data=None, limits=None)`

Return the center of mass and weighted standard deviation for data within some limits. If `data` is None, `self.data.massdat` is used. If `limits` is None, the whole range is used.

:param data: mass data to determine center of mass

:param limits: limits to restrict the calculation

:return: com, std (center of mass, weighted standard deviation)

3.56.2.4 `def unidec.UniDec.check_badness(self)`

Check for problematic variables, such as upper bounds less than lower bounds and raise warning if found.

:return:

3.56.2.5 def unidec.UniDec.convolve_peaks (self)

Convolve Peaks with Peak Shape
:return: None

3.56.2.6 def unidec.UniDec.cross_validate (self, numcrosstot = 5)

Experimental function to perform cross validation
:param numcrosstot: Number of cross validation routines to perform
:return: mean, stddtev (mean and standard deviation of mass distribution following cross validation)

3.56.2.7 def unidec.UniDec.export_config (self, f_name = None)

Export UniDec Configuration File
:param f_name: File name, Default of None will using config.confname
:return: None

3.56.2.8 def unidec.UniDec.export_params (self, e)

Export a number of different parameters about the peaks into different text files.
:param e: if e is "PostFit", it will output mass fit parameters as well
:return: None

3.56.2.9 def unidec.UniDec.fit_all_masses (self)

Fit all masses to a series of peaks, with initial guesses defined by the peak parameters.
:return: self.massfitdat, self.massfit (fit to data, fit parameters)

3.56.2.10 def unidec.UniDec.get_charge_peaks (self)

Determines total charge distribution. Imports each charge state as a peak in self.pks.
Will overwrite mass peaks.
:return: cpeaks (Z x 2 array of (charge state, intensity))

3.56.2.11 def unidec.UniDec.initialize (self)

Initialize Config, DataContainer, and Peaks
:return: None

3.56.2.12 def unidec.UniDec.integrate (self, limits, data = None)

Trapezoid ntegrate data between limits[0] and limits[1]
:param limits: [min,max] list of lower and upper bounds on integration
:param data: N x 2 array of data (mass, intensity)
If data is None (default), self.data.massdat is used.
:return: None

3.56.2.13 `def unidec.UniDec.kendrick_continuous (self, ref_mass = None, centermode = 0, nbins = 50, transformmode = 0, xaxistype = 1)`

Runs continuous Kendrick analysis on self.data.massdat
 :param ref_mass: Kendrick mass. Default is self.config.kendrickmass if it is already set and >0. Otherwise, default is oligomer mass (self.config.molig)
 :param centermode: Set range for normalization 0=(0,1),1=(-0.5,0.5). Default is 0.
 :param nbins: Set mass defect axis density. Default is 50 bins.
 :param transformmode: Set type of transformation. 0=Interpolation. 1=Integration. Default is 0.
 :param xaxistype: Set x-axis dimensions. 0=Kendrick Mass Number, 1=Mass Number * Kendrick Mass. Default is 1.
 :return: mass grid, mass defect grid, intensity grid. All with shape (len(self.data.massdat),nbins)

3.56.2.14 `def unidec.UniDec.kendrick_peaks (self, kmass = None, centermode = 1)`

Run Kendrick analysis on peaks (self.pks object)
 :param kmass: Kendrick mass. Default is prior kendrick mass if it exists and is >0. Otherwise, default is oligomer mass (self.config.molig)
 :param centermode: Set range for normalization 1=(0,1),0=(-0.5,0.5)
 :return: Array of [mass,defect] for each peak in self.pks.

3.56.2.15 `def unidec.UniDec.load_config (self, f_name)`

Import UniDec Configuration File
 :param f_name: File name
 :return: None

3.56.2.16 `def unidec.UniDec.load_default (self)`

Loads config from default location set at self.config.defaultconfig
 :return: None

3.56.2.17 `def unidec.UniDec.load_state (self, load_path)`

Load UniDec state from a zip save file.

Note: save_state is located under unidectools (ud.savestate)
 :param load_path: .zip file to load
 :return: True is successful, False if failed

3.56.2.18 `def unidec.UniDec.mass_grid_to_f_grid (self)`

Convert the mass vs charge grid to a mass vs charge offset grid.

Calculates the charge offset for each (mass,charge) point, creates a new axis of regularly spaced charge offsets (oaxis), and the interpolates a new grid of (mass, offset) from oaxis, which is output as outgrid.
 :return: oaxis, outgrid: offset axis (N) and offset grid (M x N)

3.56.2.19 `def unidec.UniDec.normalize_peaks (self)`

Noamlize everything in the peaks accoring to the type set in self.config.peaknorm
 0 = No normalization
 1 = Normalize the max value to 1
 2 = Normalize the sum to 1
 :return: None

3.56.2.20 def unidec.UniDec.open_file (self, file_name, file_directory, kwargs)

Open text or mzML file. Will create _unidecfiles directory if it does not exist.

If it finds a _conf.dat file in _unidecfiles, it will import the old configuration. Otherwise, it keeps the existing configuration but resets file names.

If silent=True is passed in **kwargs, prints are suppressed.

:param file_name: Name of file to open. May be in x y or x y z text format or in mzML format. May be tab or space delimited
 :param file_directory: Directory in which filename is located
 :return: None

3.56.2.21 def unidec.UniDec.pick_peaks (self)

Detect, Normalize, and Output Peaks

:return: None

3.56.2.22 def unidec.UniDec.process_data (self, kwargs)

Process data according to parameters in config.

Checks certain parameters to make sure the limits make sense. Will accept silent=True kwarg to suppress printing.

:return: None

3.56.2.23 def unidec.UniDec.process_mass_data (self)

Apply the same parameters used to process the data to process the mass distribution. Linearization parameters are ignored, but smoothing, baseline subtraction, normalization, and intensity threshold all apply.

:return: None

3.56.2.24 def unidec.UniDec.raw_process (self, dirname, inflag=False, binsize=1)

Processes Water's Raw files into .txt using external calls to:

self.config.rawreaderpath for MS
 self.config.cdcreaderpath for IM-MS

Default files are created with the header of the .raw file plus:

_rawdata.txt for MS
 _imraw.txt for IM-MS

:param dirname: .raw directory name
 :param inflag: If True, it will put the output .txt file inside the existing .raw directory. If False, it will put the file in the same directory that contains the .raw directory
 :param binsize: Parameter for IM-MS to specify the m/z bin size for conversion. If binsize=0, the conversion will be at full resolution (which is huge), so the default is every 1 m/z.
 :return: self.config.filename, self.config.dirname (name and location of created file)

3.56.2.25 def unidec.UniDec.reset_config (self)

Resets UniDec config to default. Should not reset paths or filenames.

:return: None

3.56.2.26 `def unidec.UniDec.run_unidec (self, silent = False, efficiency = False)`

Runs UniDec.

Checks that everything is set to go and then places external call to:
 self.config.UniDecPath for MS
 self.config.UniDecIMPath for IM-MS

If successful, calls self.unidec_imports()

If not, prints the error code.

:param silent: If True, it will suppress printing the output from UniDec

:param efficiency: Passed to self.unidec_imports()

:return: out (stdout from external UniDec call)

3.56.2.27 `def unidec.UniDec.save_default (self)`

Saves existing config in default location set at self.config.defaultconfig

:return: None

3.56.2.28 `def unidec.UniDec.unidec_imports (self, efficiency = False)`

Imports files output from the UniDec core executable into self.data.

:param efficiency: If True, it will ignore the larger files to speed up the run.


:return: None

The documentation for this class was generated from the following file:

- C:/Python/UniDec/unidec.py

3.57 GUniDec.UniDecApp Class Reference

Inheritance diagram for GUniDec.UniDecApp:



Public Member Functions

- `def __init__ (self, args, kwargs)`
- `def start (self)`
- `def on_end_session (self)`
- `def quit_application (self)`
- `def init (self, args, kwargs)`
- `def import_config`
- `def export_config`
- `def on_open (self, e)`
- `def on_open_file`
- `def on_save_state`
- `def on_load_state (self, e)`
- `def on_raw_open (self, e)`
- `def on_load_conf_file (self, e)`

- def **on_save_default** (self, e)
- def **on_load_default** (self, e)
- def **on_paste_spectrum** (self, e)
- def **on_reset** (self, e)
- def **on_dataprep_button** (self, e)
- def **on_unidec_button** (self, e)
- def **after_unidec_run** (self)
- def **on_pick_peaks** (self, e)
- def **on_plot_peaks** (self, e)
- def **on_peak_errors** (self, e)
- def **check_badness** (self)
- def **warn**
- def **on_auto** (self, e)
- def **makeplot1** (self, e)
- def **makeplot2** (self, e)
- def **makeplot3** (self, e)
- def **makeplot4** (self, a)
- def **makeplot5** (self, a)
- def **makeplot6**
- def **on_plot_composite** (self, e)
- def **make_im_plots** (self)
- def **on_plot_nativeccs** (self, e)
- def **on_replot** (self, e)
- def **make_cube_plot** (self, e)
- def **on_delete** (self, e)
- def **on_charge_states** (self, e)
- def **on_plot_offsets** (self, e)
- def **plot_integral**
- def **on_integrate**
- def **on_smash** (self)
- def **on_charge_plot** (self, e)
- def **on_batch_raw** (self, e)
- def **on_manual** (self, e)
- def **on_match** (self, e)
- def **on_mass_tools**
- def **on_peak_width_tool** (self, e)
- def **on_additional_parameters** (self, e)
- def **on_unidec_path** (self, e)
- def **on_file_name** (self, e)
- def **on_data_collector** (self, e)
- def **on_import_wizard** (self, e)
- def **on_im_tools** (self, e)
- def **on_im_extract** (self, e)
- def **on_tweet** (self, e)
- def **on_kendrick** (self, e)
- def **on_2d_grid** (self, e)
- def **on_nativez_tools** (self, e)
- def **on_export_params** (self, e)
- def **on_mass_process** (self, e)
- def **on_center_of_mass** (self, e)
- def **on_zerocharge_mass** (self, e)
- def **on_fit_masses** (self, e)
- def **on_batch**
- def **on_batch2** (self, e)
- def **on_super_batch** (self, e)
- def **on_cross_validate** (self, e)
- def **on_pdf_report** (self, e)

Public Attributes

- **eng**
- **view**
- **twittercodes**

3.57.1 Detailed Description

Main UniDec GUI Application.

Presenter contains UniDec engine at `self.eng` and main GUI window at `self.view`

3.57.2 Member Function Documentation

3.57.2.1 `def GUniDec.UniDecApp.export_config (self, file_name = None)`

Get configuration from GUI and (if `file_name` is specified) write from engine to `file_name`
:param `file_name`: Path of file to save config to
:return: None

3.57.2.2 `def GUniDec.UniDecApp.import_config (self, file_name = None)`

Import configuration from file to engine (if `file_name` is specified) and from engine to GUI.
:param `file_name`: Path of file to import
:return: None

The documentation for this class was generated from the following file:

- `C:/Python/UniDec/GUniDec.py`

3.58 unidecstructure.UniDecConfig Class Reference

Public Member Functions

- `def __init__ (self)`
- `def initialize (self)`
- `def default_colormaps (self)`
- `def config_export (self, name)`
- `def config_import (self, name)`
- `def print_config (self)`
- `def default_file_names (self)`
- `def check_badness (self)`
- `def default_high_res (self)`
- `def default_zero_charge (self)`
- `def initialize_system_paths (self)`

Public Attributes

- **infname**
- **outfname**
- **mfile**
- **manualfile**
- **confname**

- ofile
- matchfile
- peaksfile
- dirname
- filename
- extension
- imflag
- publicationmode
- discreteplot
- cmap
- peakcmap
- rawflag
- detectoreffva
- mzbins
- smooth
- subbuff
- subtype
- intthresh
- minmz
- maxmz
- numit
- zzsig
- startz
- endz
- numz
- mzsiz
- psfun
- massub
- masslb
- msig
- molig
- massbins
- adductmass
- damp
- aggressiveflag
- suppression
- isotopemode
- peakwindow
- peakthresh
- peakplotthresh
- separation
- peaknorm
- error
- mtabsig
- poolflag
- nativezub
- nativezlb
- inflate
- linflag
- integratelb
- integrateub
- mindt
- maxdt
- smoothdt
- subbufdt

- **ccslb**
- **ccsub**
- **nativeccsub**
- **nativeccslb**
- **dtsig**
- **ccsbins**
- **csig**
- **pusher**
- **zout**
- **temp**
- **pressure**
- **volt**
- **to**
- **driftlength**
- **tcal1**
- **tcal2**
- **edc**
- **gasmass**
- **twaveflag**
- **batchflag**
- **procflag**
- **runtime**
- **massdatnormtop**
- **mfileflag**
- **manualfileflag**
- **kendrickmass**
- **masslist**
- **matchlist**
- **oligomerlist**
- **manuallist**
- **zoffs**
- **gridparams**
- **cmaps**
- **badtest**
- **warning**
- **system**
- **defaultUnidecDir**
- **defaultUnidecName**
- **defaultIMName**
- **UniDecPath**
- **UniDecName**
- **UniDecIMName**
- **UniDecDir**
- **UniDecIMPath**
- **rawreaderpath**
- **cdcreaderpath**
- **defaultconfig**

3.58.1 Detailed Description

Class containing all options and configurations for UniDec GUI and Program. Contains methods to export and import config to text file for running UniDec core binaries and for storing parameters for GUI.

3.58.2 Constructor & Destructor Documentation

3.58.2.1 `def unidecstructure.UniDecConfig.__init__(self)`

Initialize Everything. Set default paths and run `self.initialize`
:return: UniDecConfig object

3.58.3 Member Function Documentation

3.58.3.1 `def unidecstructure.UniDecConfig.check_badness(self)`

Test for a few things that will crash the program:
Min is greater than Max for m/z, charge, mass, native charge, ccs, native ccs, dt
Bad IM-MS calibration values.
Peak width is zero
m/z resolution is really small.
:return: None

3.58.3.2 `def unidecstructure.UniDecConfig.config_export(self, name)`

Writes config to file give in name. Typically in format: name value.

Also exports manuallist, masslist, and oligomerlist to text files.
:param name: File name to write to.
:return: None

3.58.3.3 `def unidecstructure.UniDecConfig.config_import(self, name)`

Imports configuration from txt file. Also imports masslist, manuallist, and oligomerlist.
:param name: File to import from.
:return: None

3.58.3.4 `def unidecstructure.UniDecConfig.default_colormaps(self)`

Get default matplotlib colormaps and set names to `self.cmaps`. Ignores `CMRmap` and `CMRmap_r` because the caused problems.
:return: None

3.58.3.5 `def unidecstructure.UniDecConfig.default_file_names(self)`

Sets the default file names. For things comming into and out of the program. In theory these can be modified, but it might be risky.
:return: None

3.58.3.6 `def unidecstructure.UniDecConfig.default_high_res(self)`

Sets some defaults for high resolution spectra. Leaves other values unchanged.
:return: None

3.58.3.7 `def unidecstructure.UniDecConfig.default_zero_charge(self)`

Sets some defaults for when the zero-charge mass spectrum itself is to be deisotoped. Leaves other values unchanged.
:return: None

3.58.3.8 def unidecstructure.UniDecConfig.initialize (self)

Initialize configuration parameters but not paths. Runs self.default_colormaps
:return: None

3.58.3.9 def unidecstructure.UniDecConfig.initialize_system_paths (self)

Initialize initial paths for UniDec directories
:return: None

3.58.3.10 def unidecstructure.UniDecConfig.print_config (self)


Simple debugging command to read in the config file and print out its contents.
:return: None

The documentation for this class was generated from the following file:

- C:/Python/UniDec/unidecstructure.py

3.59 unidec_tests.UniDecTest Class Reference

Inheritance diagram for unidec_tests.UniDecTest:

**Public Member Functions**

- def **test_spectra** (self, mlist=None, ilit=None, res=1000, noise=0.0, pad=2000, window=None, plot=False, restricted=True, massbins=10, psfun=0, mzsigs=None, kwargs)
- def **simpleplotter** (self, m1, m2, e1, e2, i)
- def **plot_repeats** (self, peaks)
- def **plot_test** (self, tests, vals, testdat, testmdat, peaks)
- def **test_width** (self, kwargs)
- def **test_noise** (self, kwargs)
- def **test_zwidth** (self, kwargs)
- def **test_baseline** (self, kwargs)
- def **test_ambiguity** (self, kwargs)
- def **test_density** (self, kwargs)
- def **test_psfun** (self, kwargs)
- def **test_widths** (self, kwargs)
- def **test_sigs** (self, kwargs)
- def **test_repeat** (self, kwargs)


Additional Inherited Members

The documentation for this class was generated from the following file:

- C:/Python/UniDec/unidec_tests.py

3.60 import_wizard_grid.WizardGrid Class Reference

Inheritance diagram for import_wizard_grid.WizardGrid:



Public Member Functions

- def `__init__`
- def `sneaky_resize` (self, panel_width)
- def `showPopupMenu` (self, evt)
- def `set_labels` (self, mode)
- def `EvtDriftType` (self, evt)
- def `clear_all` (self, evt)
- def `add_dataset` (self, out)
- def `next_free_row`
- def `column_labels` (self)
- def `fill_down` (self, evt)
- def `remove_row` (self, evt)

Public Attributes

- `col_header`
- `col_conv`
- `col`
- `popupID1`

3.60.1 Detailed Description

Grid for data import wizard

3.60.2 Member Function Documentation

3.60.2.1 def import_wizard_grid.WizardGrid.fill_down (self, evt)

Try to fill down the columnn

3.60.2.2 def import_wizard_grid.WizardGrid.set_labels (self, mode)

Set initial column headers based on mode
 0 - Linear
 1 - T-wave

3.60.2.3 def import_wizard_grid.WizardGrid.showPopupMenu (self, evt)


Create and display a popup menu on right-click event

The documentation for this class was generated from the following file:

- C:/Python/UniDec/import_wizard_grid.py

3.61 datacollector.XListCtrl Class Reference

Inheritance diagram for datacollector.XListCtrl:



Public Member Functions

- def **__init__**
- def **Populate**
- def **Clear** (self)
- def **AddLine**
- def **GetList** (self)
- def **GetMaxes** (self)

Public Attributes


- **currentItem**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/datacollector.py

3.62 datacollector.YListCtrl Class Reference

Inheritance diagram for datacollector.YListCtrl:



Public Member Functions

- def **__init__**
- def **Populate**
- def **Clear** (self)

- def **AddLine**
- def **GetList** (self)

The documentation for this class was generated from the following file:

- C:/Python/UniDec/datacollector.py

3.63 nativez.zoffset Class Reference

Public Member Functions

- def **__init__** (self, args)
- def **Make** (self, offset, intensity, index, color, marker)

Public Attributes

- **offset**
- **intensity**
- **index**
- **color**
- **marker**
- **id**
- **width**
- **nstate**
- **extractwidth**

The documentation for this class was generated from the following file:

- C:/Python/UniDec/nativez.py

3.64 ZoomBox.ZoomBox Class Reference

Public Member Functions

- def [__init__](#)
- def **new_axes**
- def **update_background** (self, event)
- def **ignore** (self, event)
- def **press** (self, event)
- def **release** (self, event)
- def **update** (self)
- def **onmove** (self, event)
- def [set_active](#) (self, active)
- def [get_active](#) (self)

Public Attributes

- **crossoverpercent**
- **axes**
- **canvas**
- **visible**
- **cids**
- **active**
- **to_draw**
- **background**
- **onselect**
- **onmove_callback**
- **useblit**
- **minspanx**
- **minspany**
- **integrate**
- **smash**
- **validButtons**
- **spancoords**
- **eventpress**
- **eventrelease**
- **data_lims**
- **rectprops**
- **buttonDown**
- **prev**

3.64.1 Detailed Description

Select a min/max range of the x axes for a matplotlib Axes

Example usage::

```
from matplotlib.widgets import RectangleSelector
from pylab import *

def onselect(xmin, xmax, value, ymin, ymax):
    'eclick and erelease are matplotlib events at press and release'
    print ' x,y min position : (%f, %f)' % (xmin, ymin)
    print ' x,y max position  : (%f, %f)' % (xmax, ymax)
    print ' used button      : ', eclick.button

def toggle_selector(event):
    print ' Key pressed.'
    if event.key in ['Q', 'q'] and toggle_selector.RS.active:
        print ' RectangleSelector deactivated.'
        toggle_selector.RS.set_active(False)
    if event.key in ['A', 'a'] and not toggle_selector.RS.active:
        print ' RectangleSelector activated.'
        toggle_selector.RS.set_active(True)

x = arange(100)/(99.0)
y = sin(x)
fig = figure
axes = subplot(111)
axes.plot(x,y)

toggle_selector.RS = ZoomBox(axes, onselect, drawtype='line')
connect('key_press_event', toggle_selector)
show()
```

3.64.2 Constructor & Destructor Documentation

3.64.2.1 `def ZoomBox.ZoomBox.__init__(self, axes, onselect, drawtype = 'box', minspanx = None, minspany = None, useblit = False, lineprops = None, rectprops = None, onmove_callback = None, spancoords = 'data', button = None, data_lims = None, integrate = 0, smash = 0)`

Create a selector in axes. When a selection is made, clear the span and call onselect with

```
onselect(pos_1, pos_2)
```

and clear the drawn box/line. There pos_i are arrays of length 2 containing the x- and y-coordinate.

If minspanx is not None then events smaller than minspanx in x direction are ignored(it's the same for y).

The rect is drawn with rectprops; default
`rectprops = dict(facecolor='red', edgecolor = 'black',
alpha=0.5, fill=False)`

The line is drawn with lineprops; default
`lineprops = dict(color='black', linestyle='-',
linewidth = 2, alpha=0.5)`

Use type if you want the mouse to draw a line, a box or nothing between click and actual position by setting

```
drawtype = 'line', drawtype='box' or drawtype = 'none'.
```

spancoords is one of 'data' or 'pixels'. If 'data', minspanx and minspany will be interpreted in the same coordinates as the x and y axis, if 'pixels', they are in pixels

button is a list of integers indicating which mouse buttons should be used for rectangle selection. You can also specify a single integer if only a single button is desired. Default is None, which does not limit which button can be used.

Note, typically:

- 1 = left mouse button
- 2 = center mouse button (scroll wheel)
- 3 = right mouse button

3.64.3 Member Function Documentation

3.64.3.1 `def ZoomBox.ZoomBox.get_active(self)`

to get status of active mode (boolean variable)

3.64.3.2 `def ZoomBox.ZoomBox.set_active(self, active)`

Use this to activate / deactivate the RectangleSelector

```
from your program with an boolean variable 'active'.
```

The documentation for this class was generated from the following file:

- C:/Python/UniDec/ZoomBox.py

3.65 ZoomSpan.ZoomSpan Class Reference

Public Member Functions

- def `__init__`
- def `new_axes` (self, axes)
- def `update_background` (self, event)
- def `ignore` (self, event)
- def `press` (self, event)
- def `release` (self, event)
- def `update` (self)
- def `onmove` (self, event)

Public Attributes

- `axes`
- `canvas`
- `visible`
- `cids`
- `rect`
- `background`
- `pressv`
- `rectprops`
- `onselect`
- `onmove_callback`
- `useblit`
- `minspan`
- `buttonDown`
- `prev`
- `data_lims`

3.65.1 Detailed Description

Expansion of matplotlib embed in wx example by John Bender and Edward Abraham, see http://www.scipy.org/Matplotlib_figure_in_a_wx_panel

This version allows the user to zoom in on the figure using either a span selector or a box selector. You can also set a persistent span selector that acts as cursor references on top of whatever is plotted

ZoomSpan based on matplotlib.widgets.SpanSelector
 CursorSpan based on matplotlib.widgets.SpanSelector
 BoxZoom based on matplotlib.widgets.RectangleSelector

Brian J. Soher, Duke University, 20 October, 2010

Select a min/max range of the x or y axes for a matplotlib Axes

Example usage:

```
axes = subplot(111)
axes.plot(x,y)

def onselect(vmin, vmax):
    print vmin, vmax
span = ZoomSpan(axes, onselect, 'horizontal')
```

onmove_callback is an optional callback that will be called on mouse move with the span range

3.65.2 Constructor & Destructor Documentation

3.65.2.1 `def ZoomSpan.ZoomSpan.__init__(self, axes, onselect, minspan = None, useblit = False, rectprops = None, onmove_callback = None)`

Create a span selector in axes. When a selection is made, clear the span and call onselect with

```
onselect(vmin, vmax)
```

If minspan is not None, ignore events smaller than minspan

The span rect is drawn with rectprops; default
`rectprops = dict(facecolor='red', alpha=0.5)`

set the visible attribute to False if you want to turn off the functionality of the span selector

The documentation for this class was generated from the following file:

- C:/Python/UniDec/ZoomSpan.py

Index

- `__init__`
 - `mainwindow::Mainwindow`, 27
 - `nativez::ColorList`, 9
 - `peakstructure::Peak`, 39
 - `peakstructure::Peaks`, 40
 - `unidec::DataContainer`, 15
 - `unidec::UniDec`, 59
 - `unidecstructure::UniDecConfig`, 68
 - `ZoomBox::ZoomBox`, 74
 - `ZoomSpan::ZoomSpan`, 76
- `add_children`
 - `import_wizard_treectrl::TreeCtrlPanel`, 56
- `add_peaks`
 - `peakstructure::Peaks`, 40
- `add_root`
 - `import_wizard_treectrl::TreeCtrlPanel`, 56
- `autocorrelation`
 - `unidec::UniDec`, 59
- `autointegrate`
 - `unidec::UniDec`, 59
- `center_of_mass`
 - `unidec::UniDec`, 59
- `check_badness`
 - `unidec::UniDec`, 59
 - `unidecstructure::UniDecConfig`, 68
- `clear_all_plots`
 - `mainwindow::Mainwindow`, 27
- `ColorPlot.ColorPlot2D`, 9
- `config_export`
 - `unidecstructure::UniDecConfig`, 68
- `config_import`
 - `unidecstructure::UniDecConfig`, 68
- `convolve_peaks`
 - `unidec::UniDec`, 59
- `cross_validate`
 - `unidec::UniDec`, 60
- `datacollector.DataCollector`, 12
- `datacollector.ListCtrlPanel`, 23
- `datacollector.NetworkFrame`, 38
- `datacollector.XListCtrl`, 71
- `datacollector.YListCtrl`, 71
- `default_colormaps`
 - `unidecstructure::UniDecConfig`, 68
- `default_file_names`
 - `unidecstructure::UniDecConfig`, 68
- `default_high_res`
 - `unidecstructure::UniDecConfig`, 68
- `default_params`
 - `peakstructure::Peaks`, 40
- `default_zero_charge`
 - `unidecstructure::UniDecConfig`, 68
- `export_config`
 - `GUniDec::UniDecApp`, 65
 - `unidec::UniDec`, 60
- `export_file`
 - `import_wizard::ImportWizard`, 18
- `export_gui_to_config`
 - `mainwindow::Mainwindow`, 27
- `export_params`
 - `unidec::UniDec`, 60
- `export_then_load`
 - `import_wizard::ImportWizard`, 18
- `Extract2D.Extract2DPlot`, 15
- `fill_down`
 - `import_wizard_grid::WizardGrid`, 70
- `fit_all_masses`
 - `unidec::UniDec`, 60
- `GUniDec.Shell`, 48
- `GUniDec.UniDecApp`, 63
- `GUniDec::UniDecApp`
 - `export_config`, 65
 - `import_config`, 65
- `get_active`
 - `ZoomBox::ZoomBox`, 74
- `get_charge_peaks`
 - `unidec::UniDec`, 60
- `get_folder_path`
 - `import_wizard::ImportWizard`, 18
- `get_mass_defects`
 - `peakstructure::Peaks`, 40
- `IM_windows.IMToolExtract`, 19
- `IM_windows.IMTools`, 20
- `IM_windows.TestListCtrl4`, 51
- `IM_windows.TestListCtrlPanel4`, 54
- `import_config`
 - `GUniDec::UniDecApp`, 65
- `import_config_to_gui`
 - `mainwindow::Mainwindow`, 27
- `import_wizard.ImportWizard`, 18
- `import_wizard::ImportWizard`
 - `export_file`, 18
 - `export_then_load`, 18
 - `get_folder_path`, 18

- import_wizard_grid.WizardGrid, 70
- import_wizard_grid::WizardGrid
 - fill_down, 70
 - set_labels, 70
 - showPopupMenu, 70
- import_wizard_treectrl.TreeCtrl, 55
- import_wizard_treectrl.TreeCtrlPanel, 55
- import_wizard_treectrl::TreeCtrlPanel
 - add_children, 56
 - add_root, 56
- initialize
 - unidec::UniDec, 60
 - unidecstructure::UniDecConfig, 68
- initialize_system_paths
 - unidecstructure::UniDecConfig, 69
- integrate
 - unidec::UniDec, 60
- kendrick_continuous
 - unidec::UniDec, 60
- kendrick_peaks
 - unidec::UniDec, 61
- load_config
 - unidec::UniDec, 61
- load_default
 - unidec::UniDec, 61
- load_state
 - unidec::UniDec, 61
- mainwindow.Mainwindow, 23
- mainwindow::Mainwindow
 - __init__, 27
 - clear_all_plots, 27
 - export_gui_to_config, 27
 - import_config_to_gui, 27
 - menu_401_403, 28
 - on_about, 28
 - on_check_manual, 28
 - on_defaults, 28
 - on_exit, 28
 - on_flip_mode, 28
 - on_flip_tabbed, 28
 - on_flip_twave, 28
 - on_mass_list, 28
 - on_motion, 29
 - on_save_figure_dialog, 29
 - on_save_figure_eps, 29
 - on_save_figure_pdf, 29
 - on_save_figure_png, 29
 - on_save_figure_small, 29
 - save_all_figures, 29
 - setup_main_panel, 29
 - setup_menu, 30
 - setup_shortcuts, 30
 - setup_tool_tips, 30
 - shrink_all_figures, 30
 - shrink_figure, 30
- mass_grid_to_f_grid
 - unidec::UniDec, 61
- MassDefects.MassDefectWindow, 32
- MassFitter.MassFitter, 33
- MassModelFitter.mass, 31
- masstools.AutocorrWindow, 7
- masstools.CorrListCtrl, 10
- masstools.CorrListCtrlPanel, 10
- masstools.ManualSelection, 30
- masstools.MassSelection, 34
- masstools.TestListCtrl, 49
- masstools.TestListCtrl2, 49
- masstools.TestListCtrl3, 50
- masstools.TestListCtrlMatch, 51
- masstools.TestListCtrlPanel, 52
- masstools.TestListCtrlPanel2, 53
- masstools.TestListCtrlPanel3, 53
- masstools.TestListCtrlPanelMatch, 55
- menu_401_403
 - mainwindow::Mainwindow, 28
- miscwindows.AdditionalParameters, 5
- miscwindows.FileNameDialog, 16
- miscwindows.SaveFigureDialog, 47
- miscwindows.SingleInputDialog, 48
- mzMLImporter.mzMLImporter, 36
- nativez.ColorList, 8
- nativez.NativeZ, 36
- nativez.zoffset, 72
- nativez::ColorList
 - __init__, 9
- normalize_peaks
 - unidec::UniDec, 61
- on_about
 - mainwindow::Mainwindow, 28
- on_check_manual
 - mainwindow::Mainwindow, 28
- on_defaults
 - mainwindow::Mainwindow, 28
- on_exit
 - mainwindow::Mainwindow, 28
- on_flip_mode
 - mainwindow::Mainwindow, 28
- on_flip_tabbed
 - mainwindow::Mainwindow, 28
- on_flip_twave
 - mainwindow::Mainwindow, 28
- on_mass_list
 - mainwindow::Mainwindow, 28
- on_motion
 - mainwindow::Mainwindow, 29
- on_save_figure_dialog
 - mainwindow::Mainwindow, 29
- on_save_figure_eps
 - mainwindow::Mainwindow, 29
- on_save_figure_pdf
 - mainwindow::Mainwindow, 29
- on_save_figure_png
 - mainwindow::Mainwindow, 29

- on_save_figure_small
 - mainwindow::Mainwindow, 29
- open_file
 - unidec::UniDec, 61
- peaklistsort.TestListCtrl, 49
- peaklistsort.TestListCtrlPanel, 52
- peakstructure.Peak, 38
- peakstructure.Peaks, 39
- peakstructure::Peak
 - __init__, 39
- peakstructure::Peaks
 - __init__, 40
 - add_peaks, 40
 - default_params, 40
 - get_mass_defects, 40
 - score_peaks, 40
- peakwidthtools.PeakTools1d, 40
- peakwidthtools.PeakTools2d, 42
- pick_peaks
 - unidec::UniDec, 62
- plot1d.BarChart, 7
- plot1d.Plot1d, 44
- plot2d.Plot2d, 44
- plot2d.Plot2dMass, 45
- plot3d.CubePlot, 10
- PlotAnimations.AnimationWindow, 6
- PlottingWindow.PlottingWindow, 46
- print_config
 - unidecstructure::UniDecConfig, 69
- process_data
 - unidec::UniDec, 62
- process_mass_data
 - unidec::UniDec, 62
- raw_process
 - unidec::UniDec, 62
- reset_config
 - unidec::UniDec, 62
- run_unidec
 - unidec::UniDec, 62
- save_all_figures
 - mainwindow::Mainwindow, 29
- save_default
 - unidec::UniDec, 63
- score_peaks
 - peakstructure::Peaks, 40
- set_active
 - ZoomBox::ZoomBox, 74
- set_labels
 - import_wizard_grid::WizardGrid, 70
- setup_main_panel
 - mainwindow::Mainwindow, 29
- setup_menu
 - mainwindow::Mainwindow, 30
- setup_shortcuts
 - mainwindow::Mainwindow, 30
- setup_tool_tips
 - mainwindow::Mainwindow, 30
- showPopupMenu
 - import_wizard_grid::WizardGrid, 70
- shrink_all_figures
 - mainwindow::Mainwindow, 30
- shrink_figure
 - mainwindow::Mainwindow, 30
- twitter_interface.PinWindow, 43
- twitter_interface.TwitterWindow, 56
- UniFit.KDmodel, 21
- UniFit.kdstruct, 22
- unidec.DataContainer, 14
- unidec.UniDec, 58
- unidec::DataContainer
 - __init__, 15
- unidec::UniDec
 - __init__, 59
 - autocorrelation, 59
 - autointegrate, 59
 - center_of_mass, 59
 - check_badness, 59
 - convolve_peaks, 59
 - cross_validate, 60
 - export_config, 60
 - export_params, 60
 - fit_all_masses, 60
 - get_charge_peaks, 60
 - initialize, 60
 - integrate, 60
 - kendrick_continuous, 60
 - kendrick_peaks, 61
 - load_config, 61
 - load_default, 61
 - load_state, 61
 - mass_grid_to_f_grid, 61
 - normalize_peaks, 61
 - open_file, 61
 - pick_peaks, 62
 - process_data, 62
 - process_mass_data, 62
 - raw_process, 62
 - reset_config, 62
 - run_unidec, 62
 - save_default, 63
 - unidec_imports, 63
- unidec_imports
 - unidec::UniDec, 63
- unidec_tests.UniDecTest, 69
- unidecstructure.UniDecConfig, 65
- unidecstructure::UniDecConfig
 - __init__, 68
 - check_badness, 68
 - config_export, 68
 - config_import, 68
 - default_colormaps, 68
 - default_file_names, 68
 - default_high_res, 68

default_zero_charge, [68](#)
initialize, [68](#)
initialize_system_paths, [69](#)
print_config, [69](#)

ZoomBox.ZoomBox, [72](#)

ZoomBox::ZoomBox

 __init__, [74](#)
 get_active, [74](#)
 set_active, [74](#)

ZoomSpan.ZoomSpan, [74](#)

ZoomSpan::ZoomSpan

 __init__, [76](#)